

### 3. LC-MS/MS Complete Kits

- 3.1 ClinMass®  
Diagnostics .....106**
- 3.2 ClinMass® TDM  
Therapeutic Drug  
Monitoring .....122**
- 3.3 ClinSpot® DBS  
Dried Blood Spot Analysis ....152**

ClinMass® LC-MS/MS Complete Kits are intended for diagnosis and therapeutic drug monitoring (TDM) in clinical, toxicological and forensic laboratories. ClinSpot® LC-MS/MS Complete Kits are designed for the analysis of dried blood spot (DBS) specimen in newborn screening.

#### **ClinMass® TDM Kit System / ClinMass® TDM Kit System, on-line**

Additionally to complete kits, a novel kit system is provided for TDM. It is based on a universal ClinMass® TDM platform that can be used for the analysis of currently 155 different drugs and metabolites. The ClinMass® TDM Kit System is available both for the standard LC-MS/MS configuration as well as for on-line SPE analysis (TDM Kit System, on-line).



### 3.1.1 Alcoholism

## Ethylglucuronide and Ethylsulfate in Urine, Plasma and Serum

Ethylglucuronide (EtG) and Ethylsulfate (EtS) are known as relatively new specific parameters for the assessment of alcohol consumption, which dates back several days. The analysis is primarily performed from urine but may be extended to plasma and serum. For this purpose RECIPE offers a ClinMass® Complete Kit for analysis from urine, which can be expanded to plasma and serum samples (ClinMass® Supplement Kit, order no. MS8100). The analytical method features fast and simple sample preparation and guarantees a reliable analyte determination. The kit validation was established in compliance with the German GTFCh guidelines for quality assurance.

- Fast and simple sample preparation
- Method validation according to GTFCh guidelines



### Test Data (Urine)

Linearity:	EtG: 76.2 - 100000 µg/l EtS: 21.8 - 76356 µg/l
LOD:	EtG: 45.7 µg/l / EtS: 13.1 µg/l
LLOQ:	EtG: 76.2 µg/l / EtS: 21.8 µg/l
Intraassay precision:	EtG: 5.3 % / EtS: 4.6 %
Interassay precision:	EtG: 5.4 % / EtS: 4.9 %

### LC-MS/MS Parameters

#### LC-Parameters

HPLC pump:	isocratic, flow rates: 0.2 ml/min, 0.5 ml/min
Injection volume:	10 µl
Injection interval:	5 min
Column heater:	40 °C

#### MS-Parameters

Ion source:	ESI, negative
MS/MS-Mode:	MRM

### Sample Preparation

#### A: Urine

##### Dilution / Dotation with IS:

50 µl urine (calibrator, control, patient)	1000 µl Diluting Solution with Internal Standard
---	--

##### LC-MS/MS Analysis:

Inject 10 µl

#### B: Plasma / Serum

##### Precipitation / Dotation with IS:

100 µl sample (calibrator, control, patient)	300 µl Precipitant P with Internal Standard
---	---

mix, centrifuge


##### Dilution:

50 µl supernatant	500 µl Diluting Solution D
----------------------	-------------------------------

##### LC-MS/MS Analysis:

Inject 10 µl

Chromatogram of the ClinChek® Urine Control, Level II (order no. MS8081):



### Ordering Information

Order No.	Description	Quantity
<b>MS8000</b>	ClinMass® Complete Kit for Ethylglucuronide and Ethylsulfate in Urine for 400 assays	1 pce.
<b>Kit Components (order no. MS8000):</b>		
<b>MS8005</b>	Autosampler Washing Solution	1000 ml
<b>MS8010</b>	Mobile Phase	1000 ml
<b>MS8012</b>	Diluting Solution with Internal Standard	4 x MS8012 (4 x 100 ml)
<b>MS8713</b>	Urine Calibrator Set, lyophil., Level 0 - 6	7 x 1 x 1 ml
	Manual	1 pce.
<b>MS8100</b>	ClinMass® Supplement Kit for Ethylglucuronide and Ethylsulfate in Plasma and Serum for 100 assays	1 pce.
<b>Kit Components (order no. MS8100):</b>		
<b>MS8112</b>	Precipitant P with Internal Standard	30 ml
<b>MS8113</b>	Serum Calibrator Set, lyophil., Level 0 - 3	4 x 1 x 1 ml
<b>MS8121</b>	Diluting Solution D	50 ml
<b>Additional Components:</b>		
<b>MS8014</b>	Optimisation Mix	2 ml
<b>MS8413</b>	Urine Calibrator Set, lyophil., Level 0 - 3	4 x 1 x 1 ml
<b>Start Accessories:</b>		
<b>MS8030</b>	Analytical Column with test chromatogram	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS8080</b>	Urine Control, lyophil., Level I	10 x 1 ml
<b>MS8081</b>	Urine Control, lyophil., Level II	10 x 1 ml
<b>MS8082</b>	Urine Control, lyophil., Level III	10 x 1 ml
<b>MS8083</b>	Urine Control, lyophil., Level I, II, III	3 x 3 x 1 ml
<b>MS8180</b>	Serum Control, lyophil., Level I	10 x 1 ml
<b>MS8181</b>	Serum Control, lyophil., Level II	10 x 1 ml
<b>MS8182</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.1.2 Biogenic Amines

## Free Metanephries in Plasma

### ► manual sample preparation

Plasma free metanephries (PFM) are the most specific parameters for the diagnosis of extra-adrenal chromaffin tumors. The measurement of PFM (metanephrine (M), normetanephrine (NM) and 3-methoxytyramine (MT)) is challenging because of their polar nature, their low molecular weight and their low physiological concentration in human plasma.

This ClinMass® Complete Kit represents a reliable and standardised method for the determination of PFM. In contrast to most of the LC-MS/MS methods established in clinical routine analysis, this method requires no dilution or evaporation steps. Thereby, measurement of highly concentrated samples is guaranteed.

Alternatively to this method with manual sample preparation for standard LC-MS/MS systems, RECIPE offers a LC-MS/MS Complete Kit with online SPE analysis

- Highly sensitive and selective LC-MS/MS method
- 5-level calibrator set, controls in 3 levels
- Isotope-labelled internal standard compounds



### Test Data

Linearity (up to):	20970 ng/l (M) 19860 ng/l (NM) 5094 ng/l (MT)
Recovery:	95 - 98 %
LOD:	0.97 ng/l (M), 0.92 ng/l (NM), 0.94 ng/l (MT)
LLOQ:	2.91 ng/l (M), 2.76 ng/l (NM), 2.83 ng/l (MT)
Intraassay precision:	2.0 % (M), 4.4 % (NM), 3.7 % (MT)
Interassay precision:	3.1 % (M), 3.2 % (NM), 3.4 % (MT)

### LC-MS/MS Parameters

#### LC-Parameters

Binary HPLC pump:	Mobile Phase A / B (0.6 ml/min)
Injection volume:	20 - 50 µl
Injection interval:	4.0 min
Column heater:	25 °C

#### MS-Parameters

Ion-Source:	ESI positive
MS/MS-Mode:	MRM

### Sample Preparation

#### Extraction:

Apply (conditioned* sample preparation column):	
500 µl Sample (calibrator, control, patient)	500 µl Internal Standard IS

wash with  
Washing Solution W1, W2 and W3

#### Elution:

Elute with 2 x 100 µl Eluting Reagent E
---


collect the eluate

#### LC-MS/MS Analysis:

Inject 20 - 50 µl
-------------------

\* Conditioning with Washing Solution W1 and W2

Chromatogram of the ClinCal® Serum Calibrator, level 2 (order no. MS11013)



#### Ordering Information

Order No.	Description	Quantity
<b>MS1100</b>	<b>ClinMass® Complete Kit for Free Metanephrines in Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS1100):</b>		
<b>MS1105</b>	Autosampler Washing Solution	1000 ml
<b>MS1107</b>	Mobile Phase A	250 ml
<b>MS1108</b>	Mobile Phase B	1000 ml
<b>MS1112</b>	Internal Standard IS	50 ml
<b>MS1103</b>	Serum Calibrator Set, lyophil. (Level 0 - 4)	5 x 1 x 3 ml
<b>MS1104</b>	Optimisation Mix	2 ml
<b>MS1120</b>	Sample Preparation Columns	50 pcs.
<b>MS1121</b>	Washing Solution W1	400 ml
<b>MS1122</b>	Washing Solution W2	400 ml
<b>MS1123</b>	Washing Solution W3	200 ml
<b>MS1124</b>	Eluting Reagent E	40 ml
<b>Start Accessories:</b>		
<b>MS11030</b>	Analytical Column with test chromatogram	1 pce.
<b>Accessories:</b>		
<b>MS1161</b>	Elution Vials	100 pcs.
<b>ClinChek® Controls:</b>		
<b>MS11080</b>	Serum Control, lyophil., Level I	10 x 3 ml
<b>MS11081</b>	Serum Control, lyophil., Level II	10 x 3 ml
<b>MS11082</b>	Serum Control, lyophil., Level III	10 x 3 ml
<b>MS11083</b>	Serum Control, lyophil., Level I, II, III	3 x 3 x 3 ml

### 3.1.2 Biogenic Amines

## Free Metanephries in Plasma

### ► on-line Analysis

Plasma free metanephries (PFM) are the most specific parameters for the diagnosis of extra-adrenal chromaffin tumors. The measurement of PFM (metanephrine (M), normetanephrine (NM) and 3-methoxytyramine (MT)) is challenging because of their polar nature, their low molecular weight and their low physiological concentration in human plasma. Most LC-MS/MS methods used in clinical routine analysis involve time- and labour-consuming manual sample preparation for solid phase extraction (SPE). In contrast to this, the RECIPE ClinMass® Complete Kit represents the first commercially available CE-IVD method and provides a simple sample pretreatment with one-step protein precipitation and subsequent on-line SPE.

- Highly sensitive and selective LC-MS/MS method
- 5-level calibrator set, controls in 3 levels
- Isotope-labelled internal standard compounds



### Test Data

Linearity (up to):	5200 ng/l (M), 10000 ng/l (NM, MT)
Recovery:	83 - 116 %
LOD:	1.94 ng/l (M), 3.67 ng/l (NM), 0.94 ng/l (MT)
LLOQ:	5.82 ng/l (M), 11.0 ng/l (NM), 2.83 ng/l (MT)
Intraassay precision:	3.1 % (M), 4.5 % (NM), 3.4 % (MT)
Interassay precision:	3.2 % (M), 4.7 % (NM), 4.1 % (MT)

### LC-MS/MS Parameters


#### LC-Parameters

Pump 1:	SPE-Buffer (0.1 - 1.0 ml/min)
Pump 2:	Mobile Phase A / B (0.6 ml/min)
Injection volume:	50 µl
Injection interval:	6.0 min
Switching valve:	automatic switching valve, 6-port/3-channel
Column heater:	25 °C


#### MS-Parameters

Ion-Source:	ESI, positive
MS/MS-Mode:	MRM

### Sample Preparation



Chromatogram of the ClinCal® Serum Calibrator, level 2 (order no. MS11013)



#### Ordering Information

Order No.	Description	Quantity
<b>MS11000</b>	<b>ClinMass® Complete Kit for Free Metanephries in Plasma – on-line Analysis for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS11000):</b>		
<b>MS11005</b>	Autosampler Washing Solution	1000 ml
<b>MS11007</b>	Mobile Phase A	250 ml
<b>MS11008</b>	Mobile Phase B	2 x MS11008 (2 x 1000 ml)
<b>MS11009</b>	SPE-Buffer Load L	2 x MS11009 (2 x 1000 ml)
<b>MS11010</b>	SPE-Buffer Wash W	1000 ml
<b>MS11012</b>	Precipitant P with Internal Standard	30 ml
<b>MS11013</b>	Serum Calibrator Set, lyophil. (Level 0 - 4)	5 x 1 x 3 ml
<b>MS11020</b>	Sample Pretreatment Vials	2 x MS11020 (2 x 100 pcs.)
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS11014</b>	Optimisation Mix	2 ml
<b>Start Accessories:</b>		
<b>MS11030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS11031</b>	SPE-Column	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS11080</b>	Serum Control, lyophil., Level I	10 x 3 ml
<b>MS11081</b>	Serum Control, lyophil., Level II	10 x 3 ml
<b>MS11082</b>	Serum Control, lyophil., Level III	10 x 3 ml
<b>MS11083</b>	Serum Control, lyophil., Level I, II, III	3 x 3 x 3 ml

### 3.1.3 Metabolic Diseases

## Homocysteine in Plasma / Serum

Homocysteine is an established parameter for the risk assessment of degenerative vascular diseases which may lead to heart attack, stroke and thrombosis. The RECIPE ClinMass® Complete Kit allows a reliable determination of total homocysteine for the detection of hyperhomocysteinemia. After an easy sample preparation (reduction and precipitation step) the samples are analysed with an injection interval of 1.0 min only. Optionally, the sample preparation can be performed in 96-well plates.

- Injection interval: 1.0 min



### Test Data

Linearity:	0.1 - 106 µmol/l
LOD:	0.01 µmol/l
LLOQ:	0.1 µmol/l
Intraassay precision:	1.9 %
Interassay precision:	4.8 %

### LC-MS/MS Parameters

#### LC-Parameters

Pump:	isocratic (0.7 ml/min)
Injection volume:	1 µl
Injection interval:	1.0 min
Column heater:	room temperature

#### MS-Parameters

Ion source:	ESI, positive
MS/MS-Mode:	MRM

### Sample Preparation

#### Reduction:

50 µl sample (+ 50 µl Internal Standard IS)	50 µl Reagent A
--	-----------------

mix, incubate  
(room temperature, 5 min)

#### Precipitation:


200 µl Reagent B
------------------

1. mix, incubate (4 °C, 5 min)  
2. centrifuge

#### LC-MS/MS Analysis:

Inject 1 µl of the supernatant

Chromatogram (left) and mass traces (right) of the ClinCal® Serum Calibrator, Level 3 (order no. MS2013)



### Ordering Information

Order No.	Description	Quantity
<b>MS2000</b>	<b>ClinMass® Complete Kit for Homocysteine in Plasma / Serum for 300 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS2000):</b>		
<b>MS2005</b>	Autosampler Washing Solution	1000 ml
<b>MS2010</b>	Mobile Phase	800 ml
<b>MS2012</b>	Internal Standard IS, lyophil.	3 x MS2012 (3 x 5 ml)
<b>MS2013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
<b>MS2020</b>	Sample Preparation Vials	3 x MS2020 (3 x 100 pcs.)
<b>MS2021</b>	Reagent A	3 x MS2021 (3 x 5 ml)
<b>MS2022</b>	Reagent B	2 x MS2022 (2 x 30 ml)
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS2030</b>	Analytical Column with test chromatogram	1 pce.
<b>Accessories:</b>		
<b>FK6501</b>	Protective Sheets for 96-Well-Plates (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>FK6510</b>	96-Well-Plates (370 µl)	5 pcs.
<b>FK6511</b>	96-Well-Plates (1200 µl)	5 pcs.
<b>ClinChek® Controls:</b>		
<b>23080</b>	Plasma Control, lyophil., Level I	10 x 3 ml
<b>23081</b>	Plasma Control, lyophil., Level II	10 x 3 ml
<b>23082</b>	Plasma Control, lyophil., Level I, II	2 x 5 x 3 ml

### 3.1.3 Metabolic Diseases

## Steroids in Serum / Plasma

Steroid hormones are synthesised from cholesterol. The majority of these powerful molecules are produced by the adrenal glands, testes or ovaries. They fulfil a wide range of physiological functions after binding to specific cell receptors of various target organs. Mainly, they are involved in the carbohydrate, lipid and protein metabolism, stress response and in the regulation of the water and mineral balances, in reproductive processes and the development of secondary sex characteristics. Consequently, steroid hormones are also involved in a variety of diseases, like adrenal tumors, the Adrenogenital syndrome (CAH), the Cushing's syndrome or the Polycystic Ovarian Syndrome. The quantitative profiling of steroids is essential in diagnosis and evaluation of pathologic conditions and therapeutic responses.

RECIPE developed a fast and simple quantitative Steroid profiling LC-MS/MS method that covers the most relevant steroid hormones.

- Fast and simple sample preparation
- Injection interval of 6 minutes
- Includes second tier parameters for CAH testing



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary gradient pump:	Mobile Phase A / B (0.5 - 0.8 ml/min)
Injection volume:	10 - 40 µl
Injection interval:	6.0 min
Column heater:	40 °C

#### MS-Parameters

Ion source:	Turbo-ionspray (TIS) ESI, Positive
MS/MS-Mode:	MRM

### Sample Pretreatment

#### Precipitation:


75 µl Mixture IS/ P	50 µl sample (calibrator, control, patient)
------------------------	--

1. mix, precipitate
2. mix, centrifuge

#### LC-MS/MS Analysis:

Inject 10 - 40 µl

Chromatogram of the ClinChek® Serum Control, Level II (order no. MS12081):



#### Ordering Information

Order No.	Description	Quantity
<b>MS12000</b>	<b>ClinMass® Complete Kit for Steroids in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS12000):</b>		
<b>MS12005</b>	Autosampler Washing Solution	1000 ml
<b>MS12008</b>	Mobile Phase A	1000 ml
<b>MS12009</b>	Mobile Phase B	800 ml
<b>MS12012</b>	Internal Standard IS	10 ml
<b>MS12013</b>	Serum Calibrator Set, lyophil. (Level 0 - 6)	7 x 1 x 2 ml
<b>MS12020</b>	Sample Preparation Vials	100 pcs.
<b>MS12021</b>	Precipitant P	5 ml
	Manual	
<b>Additional components:</b>		
<b>MS12014</b>	Optimix	2 x 1 x 2 ml
<b>Start Accessories:</b>		
<b>MS12030</b>	Analytical Column with test chromatogram	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS12083</b>	Serum Control, lyophil., Level I, II, III	3 x 3 x 2 ml

### 3.1.4 Vitamin Status

## Methylmalonic Acid in Serum / Plasma / Urine, advanced

Methylmalonic acid (MMA) serves for the assessment of the vitamin B12 status. Vitamin B12 deficiency leads to an increased MMA concentration in serum, plasma and urine. This new RECIPE ClinMass® Complete Kit provides a fast MMA analysis featuring short sample preparation with only one pipetting step (serum, plasma) or direct injection of samples after dilution and addition of IS (urine). The subsequent LC-MS/MS analysis is performed within a 2 min injection interval only while the isobar component succinic acid, which may interfere with MMA, is separated reliably. Calibration and quality control is performed in a convenient and cost-efficient way. The method is fully validated using only one (serum) calibration and control set for all three sample matrices (ClinCal® 4-Level Calibrator Set and ClinChek® Controls).

- One-step sample preparation (serum / plasma)



### Test Data

Linearity:

– Serum, Plasma:	25 - 67745 nmol/l
– Urine:	651 - 1426200 nmol/l

LOD:

– Serum, Plasma:	15 nmol/l
– Urine:	391 nmol/l

LLOQ:

– Serum, Plasma:	25 nmol/l
– Urine:	651 nmol/l

Intraassay precision: 4.3 %

Interassay precision: 5.1 %

### LC-MS/MS Parameters

#### LC-Parameters

Binary gradient pump: Mobile Phase A / B (0.6 ml/min)

Injection volume: 5 µl

Injection interval: 2.0 min

Column heater: 40 °C

#### MS-Parameters

Ion source: Turbo-ionspray (TIS) ESI, negative

MS/MS-Mode: MRM

### Sample Preparation

#### Precipitation / Addition of IS:

400 µl P Precipitant (contains internal standard)	100 µl sample* (calibrator, control, patient)
--	--


mix and centrifuge

#### LC-MS/MS Analysis:

Inject 5 µl of the sample

\* Urines are diluted (50 µl urine + 1000 µl Diluting Solution D)

Chromatogram of the ClinChek® Serum Control, Level I (order no. MS5080)



#### Ordering Information

Order No.	Description	Quantity
<b>MS5100</b>	<b>ClinMass® Complete Kit, advanced, for Methylmalonic Acid in Serum / Plasma / Urine for 300 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS5100):</b>		
<b>MS5005</b>	Autosampler Washing Solution	1000 ml
<b>MS5108</b>	Mobile Phase A	1000 ml
<b>MS5109</b>	Mobile Phase B	200 ml
<b>MS5112</b>	Precipitant P with Internal Standard	3 x MS5112 (3 x 40 ml)
<b>MS5013</b>	Serum Calibrator Set, lyophil., (Level 0 - 3)	4 x 1 x 2 ml
<b>MS5020</b>	Sample Preparation Vials	3 x MS5020 (3 x 100 pcs.)
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS5114</b>	Optimisation Mix	2 ml
<b>MS5021</b>	Diluting Solution D for Urine	50 ml
<b>Start Accessories:</b>		
<b>MS5130</b>	Analytical Column with test chromatogram	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS5080</b>	Serum Control, lyophil., Level I	10 x 2 ml
<b>MS5081</b>	Serum Control, lyophil., Level II	10 x 2 ml
<b>MS5082</b>	Serum Control, lyophil., Level I, II	2 x 5 x 2 ml

### 3.1.4 Vitamin Status

## 25-OH Vitamin D2 / D3 in Plasma and Serum

### ► on-line Analysis

The ClinMass® Complete Kit with order nos. MS7000 enables a fast on-line SPE LC-MS/MS analysis of 25-OH-Vitamin D (D2, D3) with an injection interval of 3 min only. The C-3 epimer of 25-OH vitamin D3 (3-epi-25-OH vitamin D3) is separated reliably from the analyte peak. An analytical column with a run time of 4 min (order no. MS7035) is additionally available in order to obtain baseline separation.

- Fast on-line analysis  
(injection interval: 3 min)
- Separation of the C-3 epimer
- Traceable to NIST SRM 972



### Test Data

Linearity:	1.4 - 250 µg/l (D2) 2.3 - 250 µg/l (D3)
LOD:	0.4 µg/l (D2) 0.7 µg/l (D3)
LLOQ:	1.4 µg/l (D2) 2.3 µg/l (D3)
Intraassay precision:	6.8/3.9 % (D2) 4.7/3.3 % (D3)
Interassay precision:	7.6/4.3 % (D2) 6.6/4.2 % (D3)

### LC-MS/MS Parameters

#### LC-Parameters

Pump 1:	SPE-Buffer (0.1 - 5.0 ml/min)
Pump 2:	Mobile Phase (0.5 ml/min)
Injection volume:	10 - 50 µl
Injection interval:	3 min (4 min)
Switching valve:	6-port/3-channel (order no. FK1102)
Column heater:	40 °C

#### MS-Parameters

Ion source:	Heated Nebulizer (APCI), positive
MS/MS-Mode:	MRM

### Sample Pretreatment

#### Precipitation:

150 µl mixture IS/P	50 µl sample (calibrator, control, patient)
---------------------	---


1. mix, incubate (4 °C)
2. mix, centrifuge

#### on-line Analysis:

Inject 10 - 50 µl

**left:** ClinCal® Serum Calibrator, order no. MS7013, level 2 (25-OH vitamin D2 and D3) and ClinMass® Internal Standard, order no. MS7012 (d6-25-OH vitamin D3)

**right:** Human serum with 25-OH vitamin D3 and 3-epi-25-OH vitamin D3, measured with the 3-min method (top) and 4-min method (bottom)



### Ordering Information

Order No.	Description	Quantity
<b>MS7000</b>	ClinMass® Complete Kit for 25-OH Vitamin D2 / D3 in Plasma and Serum for 300 assays	1 pce.
<b>Kit Components (order no. MS7000):</b>		
<b>MS7005</b>	Autosampler Washing Solution	1000 ml
<b>MS7009</b>	SPE-Buffer	2 x MS7009 (2 x 1000 ml)
<b>MS7010</b>	Mobile Phase	700 ml
<b>MS7012</b>	Internal Standard IS	30 ml
<b>MS7013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
<b>MS7020</b>	Sample Pretreatment Vials	3 x MS7020 (3 x 100 pcs.)
<b>MS7021</b>	Precipitant P	15 ml
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS7014</b>	Optimisation Mix	2 ml
<b>MS7015</b>	Optimisation Mix 3-epi	1 ml
<b>Start Accessories:</b>		
<b>MS7030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS7035</b>	Analytical Column 3-epi with test chromatogram (for the additional analysis of 3-epi-25-OH Vitamin D3)	1 pce.
<b>MS7031</b>	SPE-Column	1 pce.
<b>Accessories:</b>		
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS7080</b>	Serum Control, lyophil., Level I	10 x 1 ml
<b>MS7081</b>	Serum Control, lyophil., Level II	10 x 1 ml
<b>MS7082</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.1.4 Vitamin Status

## 25-OH Vitamin D2 / D3 in Plasma and Serum

- ▶ automated on-line Analysis
- ▶ automated on-line Analysis with Filter-Plates

The ClinMass® Complete Kits with order nos. MS7100 and MS7200 enables fast on-line SPE LC-MS/MS analysis of 25-OH-Vitamin D (D2, D3) with an injection interval of 3 min only. In order to meet laboratory demands in high sample throughput for Vitamin D RECIPE has developed a fully automated sample preparation kit for many liquid handling devices with the option to use 96 well filter plates as well. The C-3 epimer of 25-OH vitamin D3 (3-epi-25-OH vitamin D3) is separated reliably from the analyte peak. An analytical column with a run time of 4 min (order no. MS7035) is additionally available in order to obtain baseline separation.

- Automated sample preparation
- Fast on-line analysis (injection interval: 3 min)
- Separation of the C-3 epimer
- Traceable to NIST SRM 972.



### Test Data

Linearity:	1.4 - 250 µg/l (D2) 2.3 - 250 µg/l (D3)
LOD:	0.4 µg/l (D2) 0.7 µg/l (D3)
LLOQ:	1.4 µg/l (D2) 2.3 µg/l (D3)
Intraassay precision:	6.8/3.9 % (D2) 4.7/3.3 % (D3)
Interassay precision:	7.6/4.3 % (D2) 6.6/4.2 % (D3)

### LC-MS/MS Parameters

#### LC-Parameters

Pump 1:	SPE-Buffer (0.1 - 5.0 ml/min)
Pump 2:	Mobile Phase (0.5 ml/min)
Injection volume:	10 - 50 µl
Injection interval:	3 min (4 min)
Switching valve:	6-port/3-channel (order no. FK1102)
Column heater:	40 °C

#### MS-Parameters

Ion source:	Heated Nebulizer (APCI), positive
MS/MS-Mode:	MRM

### Sample Processing

#### Precipitation:

150 µl mixture IS/P	50 µl sample (calibrator, control, patient)
---------------------	---

Chromatogram see page 119

### Ordering Information

Order No.	Description	Quantity
<b>MS7100</b>	<b>ClinMass® Complete Kit for 25-OH-Vitamin D2/D3 in Plasma and Serum for 1200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS7100):</b>		
<b>MS7005</b>	Autosampler Washing Solution	4 x MS7005 (4 x 1000 ml)
<b>MS7009</b>	SPE-Buffer	8 x MS7009 (8 x 1000 ml)
<b>MS7010</b>	Mobile Phase	4 x MS7010 (4 x 700 ml)
<b>MS7012</b>	IS Internal Standard	4 x MS7012 (4 x 30 ml)
<b>MS7013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x MS7013 (MS7013: 4 x 1 x 1 ml)
<b>FK6513</b>	96-Well-Plates (600 µl)	3 x FK6513 (3 x 5 pcs.)
<b>MS7021</b>	P Precipitant	4 x MS7021 (4 x 15 ml)
	Manual	1 pce.
<b>MS7200</b>	<b>ClinMass® Complete Kit for 25-OH-Vitamin D2/D3 in Plasma and Serum with Filter-Plates for 1200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS7200):</b>		
<b>MS7005</b>	Autosampler Washing Solution	4 x MS7005 (4 x 1000 ml)
<b>MS7009</b>	SPE-Buffer	8 x MS7009 (8 x 1000 ml)
<b>MS7010</b>	Mobile Phase	4 x MS7010 (4 x 700 ml)
<b>MS7012</b>	IS Internal Standard	4 x MS7012 (4 x 30 ml)
<b>MS7013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x MS7013 (MS7013: 4 x 1 x 1 ml)
<b>FK6514</b>	96-Well-Plates (300 µl)	3 x FK6514 (3 x 5 pcs.)
<b>FK6521</b>	96-Well-Plates (1000 µl)	3 x FK6521 (3 x 5 pcs.)
<b>MS7021</b>	P Precipitant	4 x MS7021 (4 x 15 ml)
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS7014</b>	Optimisation Mix	2 ml
<b>MS7015</b>	Optimisation Mix 3-epi	1 ml
<b>Start Accessories:</b>		
<b>MS7030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS7035</b>	Analytical Column 3-epi with test chromatogram (for the additional analysis of 3-epi-25-OH-Vitamin D3)	1 pce.
<b>MS7031</b>	SPE-Column	1 pce.
<b>Accessories:</b>		
<b>FK6501</b>	Protective sheets for 96-well-plate (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>FK6514</b>	96-Well-Plates (300 µl)	5 pcs.
<b>FK7340</b>	Sealings and sieves	4 pcs. each
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>ClinChek® Controls:</b>		
<b>MS7080</b>	Serum Control, lyophil., Level I	10 x 1 ml
<b>MS7081</b>	Serum Control, lyophil., Level II	10 x 1 ml
<b>MS7082</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.1 Immunosuppressants

## Immunosuppressants in Whole Blood, *advanced*

### ► on-line Analysis

The ClinMass® Complete Kit with order no. MS1100 enables a fast on-line SPE LC-MS/MS analysis of cyclosporine A, tacrolimus, sirolimus and everolimus with an injection interval of 2 min only.

- Fast on-line analysis (injection interval: 2 min)
- No matrix effects
- 4- and 7-level calibrator set available



### Test Data

Linearity (up to):	1750 µg/l (Cyc), 72 µg/l (Tcr, Srl, Evr)
LOD:	0.011 µg/l (Cyc), 0.004 µg/l (Tcr), 0.073 µg/l (Srl), 0.126 µg/l (Evr)
LLOQ:	0.037 µg/l (Cyc), 0.014 µg/l (Tcr), 0.244 µg/l (Srl), 0.420 µg/l (Evr)
Intraassay precision:	2.5 % (Cyc), 4.6 % (Tcr), 4.8 % (Srl), 4.9 % (Evr)
Interassay precision:	6.5 % (Cyc), 4.2 % (Tcr), 4.4 % (Srl), 4.3 % (Evr)

### LC-MS/MS Parameters

#### LC-Parameters

Pump 1:	SPE-Buffer (0.1 / 2.5 ml/min)
Pump 2:	Mobile Phase (0.5 / 1.0 ml/min)
Injection volume:	20 - 50 µl
Injection interval:	2 min
Switching valve:	6-port/3-channel (order no. FK1102)
Column heater:	60 °C

#### MS-Parameters

Ion source:	ESI, positive
MS/MS-Mode:	MRM

### Sample Pretreatment

#### Precipitation

200 µl Precipitant P	20 µl Internal Standard IS	100 µl whole blood (calibrator, control, patient)
-------------------------	-------------------------------	---

1. mix and incubate (ambient temperature)

2. mix and centrifuge


#### on-line Analysis

Inject 20 - 50 µl of the supernatant



### 3.2.1 Immunosuppressants

Chromatogram (left) and mass traces (right) of the ClinCal® Whole Blood Calibrator (order no. 9933), level 4



#### Ordering Information

Order	Description	Quantity
<b>MS1100</b>	ClinMass® Complete Kit, <i>advanced</i> , for Immunosuppressants in Whole Blood for 400 assays	1 pce.
<b>Kit Components (order no. MS1100):</b>		
<b>MS1005</b>	Autosampler Washing Solution	1000 ml
<b>MS1009</b>	SPE-Buffer	2 x MS1009 (2 x 1000 ml)
<b>MS1010</b>	Mobile Phase	800 ml
<b>MS1412</b>	Internal Standard IS, lyophil.	3 x MS1412 (3 x 3 ml)
<b>MS1020</b>	Sample Pretreatment Vials	4 x MS1020 (4 x 100 pcs.)
<b>MS1021</b>	Precipitant P	80 ml
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS1006</b>	Reagent MP	8 ml
<b>MS1014</b>	Optimisation Mix 1, lyophil.	2 ml
<b>MS1115</b>	Optimisation Mix 2, lyophil.	2 ml
<b>9028</b>	Whole Blood Calibrator, lyophil. (additional level for order nos. 9033 and 9933)	2 x 2 ml
<b>9033</b>	Whole Blood Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 2 ml
<b>9933</b>	Whole Blood Calibrator Set, lyophil. (Level 0 - 6)	7 x 1 x 2 ml
<b>Start Accessories:</b>		
<b>MS1030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS1031</b>	SPE-Column	1 pce.
<b>MS1032</b>	Guard Column Holder incl. 1 Guard Column	1 pce.
<b>MS1033</b>	Guard Column	5 pcs.
<b>Accessories:</b>		
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves for order no. FK7400	4 pcs. each
<b>ClinChek® Controls:</b>		
<b>8830</b>	Whole Blood Control, lyophil., Level I	5 x 2 ml
<b>8831</b>	Whole Blood Control, lyophil., Level II	5 x 2 ml
<b>8832</b>	Whole Blood Control, lyophil., Level III	5 x 2 ml
<b>8833</b>	Whole Blood Control, lyophil., Level I, II, III	3 x 2 x 2 ml
<b>8903</b>	Whole Blood Control, lyophil., Level IV, V	2 x 2 x 2 ml

### 3.2.1 Immunosuppressants

## Immunosuppressants in Whole Blood, *advanced*

- ▶ automated on-line Analysis
- ▶ automated on-line Analysis with Filter-Plates

The ClinMass® Complete Kits with order nos. MS1200 and MS1300 enable fast on-line SPE LC-MS/MS analysis of cyclosporine A, tacrolimus, sirolimus and everolimus with an injection interval of 2 min only. In order to meet laboratory demands in high sample throughput for Immunosuppressants Recipe has developed a fully automated sample preparation kit for many liquid handling devices with the option to use 96 well filter plates as well.

- Fully automated sample preparation
- Fast on-line analysis (injection interval: 2 min)
- No matrix effects
- 4- and 7-level calibrator set available



## Test Data

Linearity (up to):	1750 µg/l (Cyc), 72 µg/l (Tcr, Srl, Evr)
LOD:	0.011 µg/l (Cyc), 0.004 µg/l (Tcr), 0.073 µg/l (Srl), 0.126 µg/l (Evr)
LLOQ:	0.037 µg/l (Cyc), 0.014 µg/l (Tcr), 0.244 µg/l (Srl), 0.420 µg/l (Evr)
Intraassay precision:	2.5 % (Cyc), 4.6 % (Tcr), 4.8 % (Srl), 4.9 % (Evr)
Interassay precision:	6.5 % (Cyc), 4.2 % (Tcr), 4.4 % (Srl), 4.3 % (Evr)

## LC-MS/MS Parameters

LC-Parameters	
Pump 1:	SPE-Buffer (0.1 / 2.5 ml/min)
Pump 2:	Mobile Phase (0.5 / 1.0 ml/min)
Injection volume:	20 - 50 µl
Injection interval:	2 min
Switching valve:	6-port/3-channel (order no. FK1102)
Column heater:	60 °C
MS-Parameters	
Ion source:	ESI, positive
MS/MS-Mode:	MRM

## Sample Processing

### Precipitation:

150 µl mixture IS/P	50 µl sample (calibrator, control, patient)
---------------------	---



Chromatogram see page 123

### Ordering Information

Order No.	Description	Quantity
<b>MS1200</b>	<b>ClinMass® Complete Kit, advanced, for Immunosuppressants in Whole Blood for 1200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS1200):</b>		
<b>MS1005</b>	Autosampler Washing Solution	3 x MS1005 (3 x 1000ml)
<b>MS1009</b>	SPE-Buffer	6 x MS1009 (6 x 1000 ml)
<b>MS1010</b>	Mobile Phase	3 x MS1010 (3 x 800 ml)
<b>MS1412</b>	IS Internal Standard, lyophil.	9 x MS1412 (9 x 3 ml)
<b>FK6511</b>	96-Well Plates (1200 µl)	3 x FK6511 (3 x 5 pcs.)
<b>MS1021</b>	P Precipitant	3 x MS1021 (3 x 80 ml)
	Manual	1 pce.
<b>MS1300</b>	<b>ClinMass® Complete Kit, advanced, for Immunosuppressants in Whole Blood with Filter-Plates for 1200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS1300):</b>		
<b>MS1005</b>	Autosampler Washing Solution	3 x MS1005 (3 x 1000 ml)
<b>MS1009</b>	SPE-Buffer	6 x MS1009 (6 x 1000 ml)
<b>MS1010</b>	Mobile Phase	3 x MS1010 (3 x 800 ml)
<b>MS1412</b>	IS Internal Standard, lyophil.	9 x MS1412 (9 x 3 ml)
<b>FK6514</b>	96-Well-Plates (300 µl)	3 x FK6514 (3 x 5 pcs.)
<b>FK6521</b>	96-Well-Filter-Plates (1000 µl)	3 x FK6521 (3 x 5 pcs.)
<b>MS1021</b>	P Precipitant	3 x MS1021 (3 x 80 ml)
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS1014</b>	Optimisation Mix 1, lyophil.	2 ml
<b>MS1115</b>	Optimisation Mix 2, lyophil.	2 ml
<b>9028</b>	Whole Blood Calibrator, lyophil. (additional level for order nos. 9033 and 9933)	2 x 2 ml
<b>9033</b>	Whole Blood Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 2 ml
<b>9933</b>	Whole Blood Calibrator Set, lyophil. (Level 0 - 6)	7 x 1 x 2 ml
<b>Start Accessories:</b>		
<b>MS1030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS1031</b>	SPE-Column	1 pce.
<b>MS1032</b>	Guard Column Holder incl. 1 Guard Column	1 pce.
<b>MS1033</b>	Guard Columns	5 pcs.
<b>Accessories:</b>		
<b>FK6501</b>	Protective Sheets for 96-Well Plate (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>FK6514</b>	96-Well Plates (300 µl)	5 pcs.
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves for order no. FK7400	4 pcs. each
<b>ClinChek® Controls:</b>		
<b>8830</b>	Whole Blood Control, lyophil., Level I	5 x 2 ml
<b>8831</b>	Whole Blood Control, lyophil., Level II	5 x 2 ml
<b>8832</b>	Whole Blood Control, lyophil., Level III	5 x 2 ml
<b>8833</b>	Whole Blood Control, lyophil., Level I, II, III	3 x 2 x 2 ml
<b>8903</b>	Whole Blood Control, lyophil., Level IV, V	2 x 2 x 2 ml

### 3.2.2 TDM Kit System


## Antidepressants in Serum / Plasma

In addition to the classical TCAs a variety of next generation antidepressants (ADP) is available for the treatment of depressions, such as selective serotonin and norepinephrine reuptake inhibitors (SSRIs, SNRIs, NRIs) and TeCAs.

The ClinMass® Add-on Set for Antidepressants (order no. MS9400) contains all analyte specific components required for the quantification of 37 different ADPs and metabolites.

The ClinMass® TDM Platform (order no. MS9000) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (155 compounds).

- Simple and timesaving sample preparation
- Flexible and fast LC-MS/MS Analysis (3.7 min) extended to 37 substances
- IS contains 28 isotope-labeled analytes



## Test Data

Test data are available on request.

## LC-MS/MS Parameters

### LC-Parameters

Binary Gradient

Pump: Mobile Phases A/B  
(0.7 ml/min)

Injection volume: 1 - 10 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 3.7 min

Column heater: 40 °C


### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

## Sample Preparation

100 µl Precipitant P (contains IS)\*




One step!

\* Reconstitution of IS (order no. MS9412)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9413)



\*erythro- and threo-Dihydro-Bupropion can be separated by a modified gradient in an extra run

### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS9400</b>	<b>ClinMass® Add-on Set for Antidepressants in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9400):</b>		
<b>MS9412</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS9413</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9414</b>	Optimisation Mix, lyophil.	2 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9482</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.2 TDM Kit System


## Antiepileptic Drugs in Serum / Plasma

Antiepileptic Drugs (AED) are used for the treatment of epilepsy, a chronic neurological disorder, which is characterised by recurrent epileptic seizures.

The ClinMass® Add-on Set for Antiepileptic Drugs (order no. MS9200) contains all analyte specific components required for the quantification of 26 different AEDs and metabolites.

The ClinMass® TDM Platform (order no. MS9000) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System. The analytical column is suitable for the whole analyte panel (155 compounds).

- Simple and timesaving sample preparation
- Flexible and fast LC-MS/MS Analysis (4.5 min)



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary Gradient

Pump: Mobile Phases A/B  
(0.6 ml/min)

Injection volume: 1 - 10 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 4.5 min


Column heater: 40 °C

#### MS-Parameters

Ion source: ESI, positive / negative

MS/MS mode: MRM

### Sample Preparation




One step!

\* Reconstitution of IS (order no. MS9212)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9213)



### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS9200</b>	<b>ClinMass® Add-on Set for Antiepileptic Drugs in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9200):</b>		
<b>MS9212</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS9213</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9214</b>	Optimisation Mix, lyophil.	5 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9282</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml


### 3.2.2 TDM Kit System

## Antimycotics in Serum / Plasma

Invasive fungal infections, predominantly aspergillosis and candidiasis are the most important causes of morbidity and mortality in immunocompromised patients. Antimycotics are used for the treatment and prophylaxis of fungal infections. Antifungal drugs show large intra- and interindividual variabilities in pharmacokinetics, thus, the accurate measurement of the individual drug concentration in serum/plasma is essential for therapy control and success.

The ClinMass Add-On Set für Antimycotics (order no. MS9600) contains all analyte specific components required for the determination of 11 drugs and metabolites. Additionally, the determination of the active metabolite Voriconazole-N-Oxide is possible as well.

- Covers 11 antimycotics
- Includes Amphotericin B, Anidulafungin, Micafungin & Isavuconazole
- Isotopically labeled Internal Standards



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary Gradient  
Pump: Mobile Phases A/B  
(0.6 ml/min)

Injection volume: 1 – 10 µl


Injection interval: 3.5 min

Column heater: 40 °C

#### MS-Parameters

Ion source: ESI  
MS/MS mode: MRM

### Sample Preparation




**One step!**

\* Reconstitution of IS (order no. MS9412)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9613)



### Ordering Information


Order No.	Description	Quantity
MS9000	ClinMass® TDM Platform for 400 assays	1 pce.
<b>Kit Components (order no. MS9000):</b>		
MS9005	Autosampler Washing Solution	1000 ml
MS9007	Mobile Phase A	2 x MS9007 (2x 1000 ml)
MS9008	Mobile Phase B	1000 ml
MS9020	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
MS9021	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
MS9030	Analytical Column with test chromatogram	1 pce.
MS9032	Prefilter	1 pce.
<b>Accessories:</b>		
MS9022	Diluting Solution D	50 ml
MS9600	ClinMass® Add-on Set for Antimycotics in Serum / Plasma for 200 assays	1 pce.
<b>Kit Components (order no. MS9100):</b>		
MS9612	Internal Standard IS, lyophil.	5 x 5 ml
MS9613	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	
<b>Start Accessories:</b>		
MS9614	Optimisation Mix, lyophil.	2 ml
<b>ClinChek® Controls:</b>		
MS9682	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.2 TDM Kit System

## Benzodiazepines in Serum / Plasma

Benzodiazepines (BZP) are a group of psychoactive drugs with a broad range of therapeutic effects and belong to the most frequently prescribed drugs all over the world. The ClinMass® Add-on Set for Benzodiazepines (order no. MS9500) contains all analyte specific components required for the quantification of 36 different drugs and metabolites. The ClinMass® TDM Platform (order no. MS9000) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (155 compounds).

- One step sample preparation
- Flexible and fast LC-MS/MS Analysis (6.8 min) of 36 Benzodiazepines
- Z-drugs included



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary Gradient

Pump: Mobile Phases A/B  
(0.6 ml/min)

Injection volume: 1 - 10 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 6.8 min

Column heater: 40 °C


#### MS-Parameters

Ion source: ESI, positive

MS/MS-Mode: MRM

### Sample Preparation

100 µl Precipitant P (contains IS)\*




**One step!**

\* Reconstitution of IS (order no. MS9512)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS6013)



### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS9500</b>	<b>ClinMass® Add-on Set for Benzodiazepines in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9500):</b>		
<b>MS9512</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS6013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS6015 -</b>	Optimisation Mix 1 - 5, lyophil.	2 ml each
<b>MS6019</b>		
<b>ClinChek® Controls:</b>		
<b>MS6082</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.2 TDM Kit System


## Mycophenolic Acid in Serum / Plasma

Mycophenolic Acid (MPA) is widely used in immunosuppressive regimens after heart, liver and especially kidney transplantation to prevent organ rejection as well as for treatment of autoimmune diseases.

The ClinMass Add-On Set für Mycophenolic Acid (order no. MS99100) contains all analyte specific components required for the determination of the immunosuppressive drug. Additionally, RECIPE's analytical method ensures an adequate chromatographic separation from the glucuronide metabolite MPAG. Thus false positive MPA values caused by ion-source fragmentation of MPAG are avoided.

The ClinMass® TDM-Platform (order no. MS9000) contains all necessary components for sample preparation and LC-MS/MS analysis and can be used with all ClinMass® Add-On sets of the TDM Kit system in combination with the unique analytical column for all TDM panels.

- Simple and timesaving sample preparation
- Isotopically labelled Internal Standard



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary Gradient  
Pump: Mobile Phases A/B  
(0.5 ml/min)

Injection volume: 5 – 50 µl


Injection interval: 1.9 min

Column heater: 40 °C

#### MS-Parameters

Ion source: ESI  
MS/MS mode: MRM


### Sample Preparation



\* Reconstitution of IS (order no. MS9412)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Plasma Calibrator, Level 2 (order no. MS99113), spiked with MPAG



### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessory:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS99100</b>	<b>ClinMass® Add-on Set for Mycophenolic Acid in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS99100):</b>		
<b>MS99112</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS99113</b>	Serum Calibrator Set, lyophil., Level 0 - 3	4 x 1 x 1 ml
<b>MS9022</b>	Diluting Solution D	2 x MS9022 (2 x 50 ml)
<b>Start Accessories:</b>		
<b>MS99114</b>	Optimisation Mix, lyophil.	2 ml
<b>ClinChek® Controls:</b>		
<b>MS99183</b>	Plasma Control, lyophil., Level I, II, III	3 x 3 x 1 ml

### 3.2.2 TDM Kit System


## Neuroleptics in Serum / Plasma

Neuroleptics belong to the class of psychoactive drugs and are administered in the therapy of psychoses, in particular the symptomatic treatment of schizophrenia.

The ClinMass® Add-On Set for Neuroleptics (order no. MS9300) contains all analyte-specific components for all determination of 28 different neuroleptics.

The ClinMass® TDM Platform (order no. MS9000) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (155 compounds).

- Simple and timesaving sample preparation
- Flexible and fast LC-MS/MS Analysis (6.0 min)
- IS contains 25 isotope labelled analytes



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary HPLC  
pump: Mobile Phases A / B  
(0.65 ml/min)

Injection volume: 1 - 10 µl


Injection interval: 6.0 min

Column heater: 40 °C

#### MS-Parameters

Ion source: ESI  
MS/MS mode: MRM

### Sample Preparation




**One step!**

\* Reconstitution of IS (order no. MS9412)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9313)



### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS9300</b>	<b>ClinMass® Add-on Set for Neuroleptics in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9300):</b>		
<b>MS9312</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS9313</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	
<b>Start Accessories:</b>		
<b>MS9314</b>	Optimisation Mix, lyophil.	2 x 1 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9382</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.2 TDM Kit System


## Tricyclic Antidepressants in Serum/Plasma

Tricyclic antidepressants (TCA) are a group of psychoactive drugs which are mainly used for the therapy of endogene depressions, anxiety and pain management.

The ClinMass® Add-on Set for TCAs (order no. MS9100) contains all analyte specific components required for the quantification of 13 TCAs + 2 neuroleptics.

The ClinMass® TDM Platform (order no. MS9000) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System (BZP, TCA, ...).

- Simple and timesaving sample preparation
- Flexible and fast LC-MS/MS Analysis (3.8 min) of 15 substances
- IS contains 12 isotope labeled analytes



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Binary Gradient

Pump: Mobile Phases A/B  
(0.8 ml/min)

Injection volume: 1 - 10 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 3.8 min


Column heater: 40 °C

#### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

### Sample Preparation




**One step!**

\* Reconstitution of IS (order no. MS9512)  
with 5 ml Precipitant P (order no. MS9021)

## 3.2.2 TDM Kit System

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9113)



### Ordering Information

Order No.	Description	Quantity
<b>MS9000</b>	<b>ClinMass® TDM Platform for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9000):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9020</b>	Sample Preparation Vials	4 x MS9020 (4 x 100 pcs.)
<b>MS9021</b>	Precipitant P	2 x MS9021 (2 x 25 ml)
<b>Start Accessories:</b>		
<b>MS9030</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>MS9100</b>	<b>ClinMass® Add-on Set for Tricyclic Antidepressants in Serum / Plasma for 200 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9100):</b>		
<b>MS9112</b>	Internal Standard IS, lyophil.	5 x 5 ml
<b>MS9113</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9114</b>	Optimisation Mix, lyophil.	2 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9182</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Antidepressants in Serum / Plasma

In addition to the classical TCAs a variety of newer antidepressants (ADP) is available for the treatment of depressions, such as selective serotonin and norepinephrine reuptake inhibitors (SSRIs, SNRIs, NRIs) and TeCAs.

The ClinMass® Add-on Set for Antidepressants (order no. MS9450) contains all analyte specific components required for the quantification of 37 different ADPs and metabolites.

The ClinMass® TDM Platform (order no. MS9050) contains the components needed for on-line SPE LCMS/MS analysis. It can be used with all ClinMass® Add-on Sets of the “TDM Kit System, on-line”. The analytical column is suitable for the whole analyte panel (153 compounds).

- Flexible and fast LC-MS/MS Analysis (4.8 min)
- IS contains 28 isotope labeled analytes



## Test Data

Test data are available on request.

## LC-MS/MS Parameters

### LC-Parameters

Pump 1 (SPE): gradient program (0.1 - 3.5 ml/min)

Pump 2

(Mobile Phase A, B): gradient program (0.7 ml/min)

Injection volume: 3 - 15 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 4.8 min

Column heater: 40 °C

### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM


## on-line SPE

Direct injection of sample and Internal Standard IS in the ratio 1:2 (autosampler).

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9413)



\*erythro- and threo-Dihydro-Bupropion can be separated by a modified gradient in an extra run

#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2 x MS9059 (2 x 2500 ml)
<b>MS9060</b>	SPE-Buffer Wash W	2 x MS9060 (2 x 2500 ml)
<b>MS9071</b>	Reconstitution Solution R	5 ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves	4 pcs. each
<b>MS9450</b>	<b>ClinMass® Add-on Set, on-line, for Antidepressants in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9450):</b>		
<b>MS9462</b>	Internal Standard IS, lyophil.	5 ml
<b>MS9413</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9414</b>	Optimisation Mix, lyophil.	2 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9482</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Antiepileptic Drugs in Serum / Plasma

Antiepileptic drugs (AED) are used for the treatment of epilepsy, a chronic neurological disorder, which is characterised by recurrent epileptic seizures.

The ClinMass® Add-on Set for Antiepileptic Drugs (order no. MS9250) contains all analyte specific components required for the quantification of 27 different AEDs and metabolites.

The ClinMass® TDM Platform (order no. MS9050) contains the components needed for on-line SPE LCMS/MS analysis. It can be used with all ClinMass® Add-on Sets of the “TDM Kit System, on-line” (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (153 compounds).

- No manual sample preparation needed
- Flexible and fast LC-MS/MS Analysis (7.0 min)
- IS contains 18 isotope labeled analytes



## Test Data

Test data are available on request.

## LC-MS/MS Parameters

### LC-Parameters

Pump 1 (SPE): gradient program (0.1 - 4 ml/min)

Pump 2

(Mobile Phase A, B): gradient program (0.6 ml/min)

Injection volume: 3 - 15 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 7.0 min

Column heater: 40 °C

### MS-Parameters

Ion source: ESI, positive / negative

MS/MS mode: MRM

## Sample Preparation


No manual sample preparation needed!

Direct injection of sample and Internal Standard IS in the ratio 1:2 (autosampler).

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9213)



#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2 x MS9059 (2 x 2500 ml)
<b>MS9060</b>	SPE-Buffer Wash W	2 x MS9060 (2 x 2500 ml)
<b>MS9071</b>	Reconstitution Solution R	5 ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves	4 pcs. each
<b>MS9250</b>	<b>ClinMass® Add-on Set, on-line, for Antiepileptic Drugs in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9250):</b>		
<b>MS9262</b>	Internal Standard IS, lyophil.	5 ml
<b>MS9213</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9214</b>	Optimisation Mix, lyophil.	2 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9282</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Antimycotics in Serum / Plasma

Invasive fungal infections, predominantly aspergillosis and candidiasis are the most important causes of morbidity and mortality in immunocompromised patients. Antimycotics are used for the treatment and prophylaxis of fungal infections. Antifungal drugs show large intra- and interindividual variabilities in pharmacokinetics, thus, the accurate measurement of the individual drug concentration in serum/plasma is essential for therapy control and success.

The ClinMass Add-On Set für Antimycotics (order no. MS9600) contains all analyte specific components required for the determination of 10 drugs and metabolites. Additionally, the determination of the active metabolite Voriconazole-N-Oxide is possible as well.

The ClinMass® TDM-Platform (order no. MS9000) contains all necessary components for sample preparation and LC-MS/MS analysis and can be used with all ClinMass® Add-On sets of the TDM Kit system in combination with the unique analytical column for all TDM panels.

- No manual sample preparation needed
- Flexible and fast LC-MS/MS Analysis: 3.5 min
- IS contains 7 isotope labelled analytes



## Test Data

Test data are available on request.

## LC-MS/MS Parameters

### LC-Parameters

Pump 1 (SPE): gradient program (2.0 -4.0 ml/min)

Pump 2

(Mobile Phase A, B): gradient program (0,05 ml/min)

Injection volume: 3-15 µl\*

Injection interval: 3.5 min

Column heater: 40 °C

\* depending on the sensitivity of the mass spectrometer

### MS-Parameters

Ion source: ESI, positive/negative

MS/MS mode: MRM

## Sample Preparation


No manual sample preparation needed!

Direct injection of sample and Internal Standard IS in the ratio 1:1 (autosampler).

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9613)



#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. 9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1 x 1000 ml
<b>MS9007</b>	Mobile Phase A	2 x 1000 ml
<b>MS9008</b>	Mobile Phase B	1 x 1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2 x 2500 ml
<b>MS9060</b>	SPE-Buffer Wash W	2 x 2500 ml
<b>MS9071</b>	Reconstitution Solution R	1 x 5ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves for order no. FK7400	4 pcs. each
<b>MS9650</b>	<b>ClinMass® Add-on Set, on-line, for Antimycotics in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9650):</b>		
<b>MS9662</b>	Internal Standard IS, lyophil.	5 ml
<b>MS9613</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4x1x1 ml
	Manual	1 pce
<b>Separately available components:</b>		
<b>MS9662</b>	Internal Standard IS, lyophil.	5 ml
<b>MS9613</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
<b>Start Accessories:</b>		
<b>MS9614</b>	Optimisation Mix, lyophil.	2 ml
<b>ClinChek® Controls:</b>		
<b>MS9682</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Benzodiazepines in Serum / Plasma

Benzodiazepines are a group of psychoactive drugs with a broad range of therapeutic effects and belong to the most frequently prescribed drugs all over the world.

The ClinMass® Add-on Set for Benzodiazepines (order no. MS9550) contains all analyte specific components required for the quantification of 36 different drugs and metabolites.

The ClinMass® TDM Platform (order no. MS9050) contains the components needed for on-line SPE LCMS/MS analysis. It can be used with all ClinMass® Add-on Sets of the “TDM Kit System, on-line”. The analytical column is suitable for the whole analyte panel (153 compounds).

- No manual sample preparation needed
- Flexible and fast LC-MS/MS Analysis (7.5 min)
- IS contains 20 isotope labeled analytes



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Pump 1 (SPE): gradient program (0.1 - 4 ml/min)

Pump 2

(Mobile Phase A, B): gradient program (0.6 ml/min)

Injection volume: 2 - 20 µl\*

Injection interval: 7.5 min

Column heater: 40 °C

\* depending on the sensitivity of the mass spectrometer

#### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

### Sample Preparation


No manual sample preparation needed!

Direct injection of sample and Internal Standard IS in the ratio 1:1 (autosampler).

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS6013)



#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2 x MS9059 (2 x 2500 ml)
<b>MS9060</b>	SPE-Buffer Wash W	2 x MS9060 (2 x 2500 ml)
<b>MS9071</b>	Reconstitution Solution R	5 ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves	4 pcs. each
<b>MS9550</b>	<b>ClinMass® Add-on Set, on-line, for Benzodiazepines in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9550):</b>		
<b>MS9562</b>	Internal Standard IS, lyophil.	5 ml
<b>MS6013</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS6015 -</b>	Optimisation Mix 1 - 5, lyophil.	2 ml each
<b>MS6019</b>		
<b>ClinChek® Controls:</b>		
<b>MS6082</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Neuroleptics in Serum / Plasma

Neuroleptics belong to the class of psychoactive drugs and are administered in the therapy of psychoses, in particular the symptomatic treatment of schizophrenia.

The ClinMass® Add-On Set, for Neuroleptics, on-line, (order no. MS9350) contains all analyte-specific components for the determination of 28 different neuroleptics. The analysis is performed by HPLC with tandem mass spectrometry (LC-MS/MS).

The ClinMass® TDM Platform (order no. MS9050) contains the components needed for sample preparation and LC-MS/MS analysis. It can be used with all ClinMass® Add-on Sets of the TDM Kit System (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (153 compounds).

- Flexible and fast on-line SPE Analysis (7.0 min)
- Comprehensive analyte spectrum
- IS contains 25 isotope labelled analytes



## Test Data

Test data are available on request.

## LC-MS/MS Parameters

### LC-Parameters

Binary HPLC pump: Mobile Phases A / B (0.65 ml/min)

Binary HPLC pump: SPE-Buffer L / W (0.1 - 4 ml/min)

Injection volume: 3 - 15 µl

Injection interval: 7.0 min

Column heater: 40 °C

### MS-Parameters

Ion source: ESI

MS/MS mode: MRM

## Sample Preparation


No manual sample preparation needed!

Direct injection of sample and Internal Standard IS in the ratio 1:2 (autosampler)

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9313)



#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2500 ml
<b>MS9060</b>	SPE-Buffer Wash W	2500 ml
<b>MS9071</b>	Reconstitution Solution R	5 ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves for order no. FK7400	4 pcs. each
<b>MS9350</b>	<b>ClinMass® Add-on Set, on-line, for Neuroleptics in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9350):</b>		
<b>MS9362</b>	IS Interner Standard, lyophil.	5 ml
<b>MS9313</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	
<b>Start Accessories:</b>		
<b>MS9314</b>	Optimisation Mix, lyophil.	2 x 1 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9382</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.2.3 TDM Kit System, on-line


## Tricyclic Antidepressants in Serum / Plasma

Tricyclic antidepressants (TCA) are a group of psychoactive drugs which are mainly used for the therapy of endogene depressions, anxiety and pain management.

The ClinMass® Add-on Set for TCAs (order no. MS9150) contains all analyte specific components required for the quantification of 13 TCAs and 2 atypical neuroleptics.

The ClinMass® TDM Platform (order no. MS9050) contains the components needed for on-line SPE LCMS/ MS analysis. It can be used with all ClinMass® Add-on Sets of the “TDM Kit System, on-line” (BZP, TCA, ...). The analytical column is suitable for the whole analyte panel (153 compounds).

- No manual sample preparation needed
- Flexible and fast LC-MS/MS Analysis (3.8 min)
- IS contains 12 isotope labeled analytes



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Pump 1 (SPE): gradient program (0.1 - 5 ml/min)

Pump 2

(Mobile Phase A, B): gradient program (0.8 ml/min)

Injection volume: 3 - 15 µl\*

\* depending on the sensitivity of the mass spectrometer

Injection interval: 3.8 min

Column heater: 40 °C

#### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

### Sample Preparation


No manual sample preparation needed!

Direct injection of sample and Internal Standard IS in the ratio 1:2 (autosampler).

Direct injection!

### 3.2.3 TDM Kit System, on-line

Chromatogram of the ClinCal® Serum Calibrator, Level 2 (order no. MS9113)



#### Ordering Information

Order No.	Description	Quantity
<b>MS9050</b>	<b>ClinMass® TDM Platform, on-line for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9050):</b>		
<b>MS9005</b>	Autosampler Washing Solution	1000 ml
<b>MS9007</b>	Mobile Phase A	2 x MS9007 (2 x 1000 ml)
<b>MS9008</b>	Mobile Phase B	1000 ml
<b>MS9059</b>	SPE-Buffer Load L	2 x MS9059 (2 x 2500 ml)
<b>MS9060</b>	SPE-Buffer Wash W	2 x MS9060 (2 x 2500 ml)
<b>MS9071</b>	Reconstitution Solution R	5 ml
<b>Start Accessories:</b>		
<b>MS9035</b>	Analytical Column with test chromatogram	1 pce.
<b>MS9031</b>	SPE-Column	1 pce.
<b>MS9032</b>	Prefilter	1 pce.
<b>Accessories:</b>		
<b>MS9022</b>	Diluting Solution D	50 ml
<b>FK7400</b>	Inline-Filter (stainless steel sieve, free of dead volume)	1 pce.
<b>FK7340</b>	Sealings and sieves	4 pcs. each
<b>MS9150</b>	<b>ClinMass® Add-on Set, on-line, for Tricyclic Antidepressants in Serum / Plasma for 400 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS9150):</b>		
<b>MS9162</b>	Internal Standard IS, lyophil.	5 ml
<b>MS9113</b>	Serum Calibrator Set, lyophil. (Level 0 - 3)	4 x 1 x 1 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS9114</b>	Optimisation Mix, lyophil.	2 x 2 ml
<b>ClinChek® Controls:</b>		
<b>MS9182</b>	Serum Control, lyophil., Level I, II	2 x 5 x 1 ml

### 3.3.1 Newborn Screening (NBS)

## Amino Acids and Acylcarnitines in Dried Blood Spots (DBS)

The determination of the amino acids and acylcarnitines profile is an integral part of national newborn screening programmes. For this purpose RECIPE offers an approved ClinSpot® LC-MS/MS Complete Kit for the use in accredited screening laboratories (order no. MS10000). The kit contains all reagents needed for the analysis of 26 analytes. The determination of additional acylcarnitines is possible. The entire sample preparation is performed with 96-well-plates, in order to ensure a high sample throughput.

- Method according to CLSI NBS-04

### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Pump: isocratic (0.1 ml/min)

Injection volume: 5 - 20 µl

Injection interval: 1.0 - 1.5 min

#### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

### Sample Preparation

#### Extraction:

well-plate 1: 3.2 mm disc

100 µl Internal Standard IS

extract

#### Transfer:

Transfer supernatant into well-plate 2

evaporate

#### Derivatisation:

Add 50 µl Reagent B

1. incubate (60 °C)
2. evaporate

#### Reconstitution:

Add 100 µl Reagent C

reconstitute

#### LC-MS/MS Analysis:

Inject 10 - 20 µl



### 3.3.1 Newborn Screening (NBS)

Amino acids	Acylcarnitines	Additional acylcarnitines	
Alanine	Carnitine	C3DC-Carnitine	C16:1-Carnitine
Arginine	C2-Carnitine	C4OH-Carnitine	C16:1OH-Carnitine
Aspartic acid	C3-Carnitine	C4DC-Carnitine	C16OH-Carnitine
Citrulline	C4-Carnitine	C5:1-Carnitine	C18:2-Carnitine
Glutamic acid	C5-Carnitine	C5OH-Carnitine	C18:1-Carnitine
Glycine	C5DC-Carnitine	C6DC-Carnitine	C18:2OH-Carnitine
Leucine	C6-Carnitine	C8:1-Carnitine	C18:1OH-Carnitine
Methionine	C8-Carnitine	C10:2-Carnitine	C18OH-Carnitine
Ornithine	C10-Carnitine	C10:1-Carnitine	
Phenylalanine	C12-Carnitine	C12:1-Carnitine	
Proline	C14-Carnitine	C14:2-Carnitine	
Tyrosine	C16-Carnitine	C14:1-Carnitine	
Valine	C18-Carnitine	C14OH-Carnitine	

#### Ordering Information


Order No.	Description	Quantity
<b>MS10000</b>	<b>ClinSpot® Complete Kit for Amino Acids and Acylcarnitines in Dried Blood Spots (DBS) for 960 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS10000):</b>		
<b>MS10005</b>	Autosampler Washing Solution	1000 ml
<b>MS10010</b>	Mobile Phase	1000 ml
<b>MS10012A</b>	Internal Standard IS, lyophil.	100 ml
<b>MS10021</b>	Reagent A	100 ml
<b>MS10022</b>	Reagent B	50 ml
<b>MS10023</b>	Reagent C	100 ml
<b>MS10040</b>	96-Well-Plates (370 µl)	6 x MS10040 (6 x 5 pcs.)
<b>MS10041</b>	Covers for 96-well-plates	2 x MS10041 (2 x 15 pcs.)
<b>MS10182</b>	Dried Blood Spot Control, Level I, II	2 x 1 x 3 spots
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS10012</b>	Internal Standard IS, lyophil.	25 ml
<b>MS10014</b>	Optimisation Mix	2 x 1 x 5 ml
<b>Accessories:</b>		
<b>MS10042</b>	Protective sheets for 96-well-plate (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>MS10045</b>	Backpressure regulator, PEEK version, 34 bar	1 pce.
<b>MS10046</b>	Replacement cartridge gold coat, 34 bar	1 pce.
<b>FK7400</b>	Inline-Filter	1 pce.
<b>FK7340</b>	Sealings and Sieves for order no. FK7400	4 pcs. each
<b>ClinChek® Controls:</b>		
<b>MS10182</b>	Dried Blood Spot Control, Level I, II	2 x 1 x 3 spots

### 3.3.1 Newborn Screening (NBS)

## Amino Acids and Acylcarnitines in Dried Blood Spots (DBS) with Filter-Plates

The determination of the amino acids and acylcarnitines profile is an integral part of national newborn screening programmes. For this purpose RECIPE offers an approved ClinSpot® LC-MS/MS Complete Kit for the use in accredited screening laboratories (order no. MS10100). The kit contains all reagents needed for the analysis of 26 analytes. The determination of additional acylcarnitines is possible. With this kit (MS10100) the sample preparation is performed with 96-well-plates combined with filter-plates.

- Method according to CLSI NBS-04
- Timesaving in sample preparation



### Test Data

Test data are available on request.

### LC-MS/MS Parameters

#### LC-Parameters

Pump: isocratic (0.1 ml/min)

Injection volume: 5 - 20 µl

Injection interval: 1.0 - 1.5 min

#### MS-Parameters

Ion source: ESI, positive

MS/MS mode: MRM

### Sample Preparation

#### Extraction:

filter-plate + well-plate: 3.2 mm disc	100 µl Internal Standard IS
---	--------------------------------

1. extract
2. centrifuge
3. evaporate filtrate

#### Derivatisation:

Add 50 µl Reagent B
---------------------

1. incubate (60 °C)
2. evaporate

#### Reconstitution:

Add 100 µl Reagent C
----------------------

reconstitute

#### LC-MS/MS Analysis:

Inject 10 - 20 µl
-------------------

### 3.3.1 Newborn Screening (NBS)

Amino acids	Acylcarnitines	Additional acylcarnitines	
Alanine	Carnitine	C3DC-Carnitine	C16:1-Carnitine
Arginine	C2-Carnitine	C4OH-Carnitine	C16:1OH-Carnitine
Aspartic acid	C3-Carnitine	C4DC-Carnitine	C16OH-Carnitine
Citrulline	C4-Carnitine	C5:1-Carnitine	C18:2-Carnitine
Glutamic acid	C5-Carnitine	C5OH-Carnitine	C18:1-Carnitine
Glycine	C5DC-Carnitine	C6DC-Carnitine	C18:2OH-Carnitine
Leucine	C6-Carnitine	C8:1-Carnitine	C18:1OH-Carnitine
Methionine	C8-Carnitine	C10:2-Carnitine	C18OH-Carnitine
Ornithine	C10-Carnitine	C10:1-Carnitine	
Phenylalanine	C12-Carnitine	C12:1-Carnitine	
Proline	C14-Carnitine	C14:2-Carnitine	
Tyrosine	C16-Carnitine	C14:1-Carnitine	
Valine	C18-Carnitine	C14OH-Carnitine	

#### Ordering Information

Order No.	Description	Quantity
<b>MS10100</b>	<b>ClinSpot® Complete Kit for Amino Acids and Acylcarnitines in Dried Blood Spots (DBS) with Filter-Plates, for 960 assays</b>	1 pce.
<b>Kit Components (order no. MS10100):</b>		
<b>MS10005</b>	Autosampler Washing Solution	1000 ml
<b>MS10010</b>	Mobile Phase	1000 ml
<b>MS10012A</b>	Internal Standard IS, lyophil.	100 ml
<b>MS10021</b>	Reagent A	100 ml
<b>MS10022</b>	Reagent B	50 ml
<b>MS10023</b>	Reagent C	100 ml
<b>MS10040</b>	96-Well-Plates (370 µl)	3 x MS10040 (3 x 5 pcs.)
<b>MS10041</b>	Covers for 96-well-plates	15 pcs.
<b>MS10140</b>	96-Well-Filter-Plates (500 µl) with Covers	3 x MS10140 (3 x 5 pcs.)
<b>MS10182</b>	Dried Blood Spot Control, Level I, II	2 x 1 x 3 spots
	Manual	1 pce.
<b>Additional Components:</b>		
<b>MS10012</b>	Internal Standard IS, lyophil.	25 ml
<b>MS10014</b>	Optimisation Mix	2 x 1 x 5 ml
<b>Accessories:</b>		
<b>MS10042</b>	Protective sheets for 96-well-plate (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>MS10045</b>	Backpressure regulator, PEEK version, 34 bar	1 pce.
<b>MS10046</b>	Replacement cartridge gold coat, 34 bar	1 pce.
<b>FK7400</b>	Inline-Filter	1 pce.
<b>FK7340</b>	Sealings and Sieves for order no. FK7400	4 pcs. each
<b>ClinChek® Controls:</b>		
<b>MS10182</b>	Dried Blood Spot Control, Level I, II	2 x 1 x 3 spots

### 3.3.1 Newborn Screening (NBS)

## Homocysteine in Dried Blood Spots (DBS)

The determination of homocysteine in dried blood spots (DBS) is used for the newborn screening for inherited defects of the methionine metabolism (homocystinuria). The RECIPE ClinSpot® Complete Kit allows a reliable determination of total homocysteine from DBS. After sample preparation (reduction and precipitation step) the sample analysis is performed with an injection interval of 1.0 min only. Optionally, the sample preparation can be performed in 96-well-plates.

- Injection interval: 1.0 min



### Test Data

Linearity:	0.2 - 104 µmol/l
LOD:	0.07 µmol/l
LLOQ:	0.2 µmol/l
Intraassay precision:	2.8 %
Interassay precision:	3.8 %

### LC-MS/MS Parameters

#### LC-Parameters

Pump:	isocratic (0.7 ml/min)
Injection volume:	1 µl
Injection interval:	1.0 min
Column heater:	room temperature

#### MS-Parameters

Ion source:	ESI, positive
MS/MS-Mode:	MRM

### Sample Preparation

#### Reduction:

Filter Spot Sample (+ 15 µl Internal Standard)	15 µl Reagent A
---	-----------------

mix

#### Precipitation

100 µl Reagent B
------------------


1. mix  
2. centrifuge

#### LC-MS/MS Analysis:

Inject 1 µl of the supernatant

### 3.3.1 Newborn Screening (NBS)

Chromatogram (left) and mass traces (right) of the ClinCal® Dried Blood Spot Calibrator, Level 3 (order no. MS2113)



#### Ordering Information

Order No.	Description	Quantity
<b>MS2100</b>	<b>ClinSpot® Complete Kit for Homocysteine in Dried Blood Spots (DBS) for 300 assays</b>	<b>1 pce.</b>
<b>Kit Components (order no. MS2100):</b>		
<b>MS2005</b>	Autosampler Washing Solution	1000 ml
<b>MS2010</b>	Mobile Phase	800 ml
<b>MS2112</b>	Internal Standard IS, lyophil.	5 ml
<b>MS2113</b>	Dried Blood Spot Calibrator Set (Level 0 - 3)	4 x 1 x 3 spots
<b>MS2020</b>	Sample Preparation Vials	3 x MS2020 (3 x 100 pcs.)
<b>MS2021</b>	Reagent A	5 ml
<b>MS2022</b>	Reagent B	30 ml
	Manual	1 pce.
<b>Start Accessories:</b>		
<b>MS2030</b>	Analytical Column with test chromatogram	1 pce.
<b>Accessories:</b>		
<b>FK6501</b>	Protective Sheets for 96-Well-Plates (PE/PP foil, 80 x 140 mm)	50 pcs.
<b>FK6510</b>	96-Well-Plates (370 µl)	5 pcs.
<b>ClinChek® Controls:</b>		
<b>MS2182</b>	Dried Blood Spot Control, Level I, II	2 x 1 x 3 spots