



KB03039

Cysteine Assay Kit

**96 well plate
100/200/400 tests**

Table of contents

1. General information	1
2. Technical specifications	2
3. Materials and storage	3
4. Introduction	4
5. Assay Principle	4
6. Assay preparation	5
7. Sample preparation	7
8. Assay protocol	8
9. Data analysis	9
10. Troubleshooting	11
11. Additional information	13
12. Related products	13
13. Warranties and limitation of liability	14

1. General information

PRECAUTIONS

Please read this manual carefully before beginning the assay.

This product is designed for **research use only**. It is not approved for human or animal use or clinical diagnosis. All chemicals should be handled with care and in accordance with laboratory safety practices. It is recommended to use basic Personal Protective Equipment.

Do not use after the expiration date stated on the packaging.

Do not mix or substitute reagents or materials from other kit batches or vendors.

For the **material safety data sheet** (MSDS) please contact us at info@bioquochem.com

TECHNICAL RECOMMENDATIONS

Store reagents as indicated in **Materials and storage** section.

Be sure to keep the bottle capped when not in use.

Let the components reach room temperature (RT) before use.

Immediately before use, gently invert and rotate reagent bottles several times to mix the contents thoroughly.

Avoid foaming or bubbles when mixing or reconstituting components.

Avoid cross contamination of samples or reagents by changing pipette tips between sample, standard and reagent additions.

Be sure to use the optimal microplate for the assay. Flat bottom transparent microplates for UV/VIS applications, and black microplates for fluorescence measurements.

2. Technical specifications

Available sizes

100/200/400 tests

Required sample volume

20 µL/test

Compatible samples

Animal and Plant tissue homogenates, cells, bacteria, plasma, serum, urine and other biological samples

Type of detection

Colorimetric (600 nm)

3. Materials and storage

MATERIALS SUPPLIED

Store kit components as indicated below:

Item	No. Tests	Units	Storage
Reagent A	100	1	4 °C
	200	2	
	400	4	
Reagent B	100	1	4 °C
	200	2	
	400	4	
Reagent C	100	1	RT
	200	2	
	400	4	
Reagent D	100	1	4 °C
	200	1	
	400	1	
Standard	100	1	4 °C
	200	1	
	400	1	
Transparent 96-Well Microplate	100	1	RT
	200	2	
	400	4	

MATERIALS NEEDED BUT NOT SUPPLIED

- Double distilled water (ddH₂O) as Milli-Q Ultrapure Water
- Incubator
- Labware materials (micropipettes, tubes, stirring/mixing equipment)
- Colorimetric microplate reader – equipped with filter for OD 600 nm

STORAGE CONDITIONS

On receipt, store kit components as indicated above. Under these conditions, the reagents are stable in the original packaging until the expiration date indicated on the outside of the box. **Reagents B and Standard are light sensitive** and should be stored in the dark. Standard solutions are unstable in the presence of oxygen. Prepare a fresh set of standards for every use.

4. Introduction

Cysteine (Cys) is the only sulfur-containing amino acid with a sulfhydryl group, an important structural and functional amino acid of proteins, which is synthesized from homocysteine in animals. Cys participates in the formation of protein disulfide bonds and it is often a component of the active center of enzymes. It also provides sulfhydryl groups for other biological reactions.

Cys can be accumulated in large amounts on the skin and mucous membranes to maintain normal metabolism of the skin during keratin and melanin production. Cys is also involved in several redox pathways and it is a precursor of glutathione. Low levels of Cys have been related to damage from reactive oxygen species (ROS), while high levels have been related to several cardiovascular and metabolic diseases.

BQC Cysteine Assay Kit is a quick, easy, and reproducible assay to quantify cysteine content in a wide variety of samples.

5. Assay Principle

This assay is based on the reduction of phosphotungstic acid by Cys to produce tungsten blue. The concentration of this blue compound can be measured at 600 nm and is proportional to the Cys content in the sample.



Principle of Cysteine Assay Kit

6. Assay preparation

REAGENT PREPARATION

All assay reagents not listed below are ready to use as supplied. Allow the reagents to reach room temperature before use.

R.B Working Solution: For 100 tests, add 2.5 mL deionized water to **Reagent B** for full dissolution, and then add 0.625 mL of **Reagent C**. Mix well, incubate in a boiling water bath for 2 h and cover to prevent water loss. After cooling, add 10 mL deionized water and equilibrate at room temperature before use.


⚠ **CAUTION:** The remaining R.B. Working Solution can be stored at 4°C for 2 weeks. Keep in the dark.

Standard Solution (10 mM): add 8.264 mL of deionized water to the **Standard** to prepare **10 mM Standard Solution**.

⚠ **CAUTION:** Standard Solution must be freshly prepared and used immediately. Keep in the dark.

STANDARD CALIBRATION

Prepare Cysteine standards for the calibration curve from the Standard Solution 10 mM making eight serial dilutions according to the following Figure and Table. Prepare the standards immediately prior to each assay. Vortex tubes thoroughly. Discard standard solutions after use.



Preparation of standards from the Standard Solution making serial dilutions

Standard	Standard Solution (μL)	ddH ₂ O (μL)	Standard Concentration (mM)
Std 1	200 Standard Solution 10 mM	0	10
Std 2	100 Std1	100	5
Std 3	100 Std2	100	2.5
Std 4	100 Std3	100	1.25
Std 5	100 Std4	100	0.625
Std 6	100 Std5	100	0.3125
Std 7	100 Std6	100	0.156
Std 8 (Reagent Blank)	0	200	0

PLATE SET UP

BQC recommends running the standards and samples at least in duplicate (triplicate recommended). There is no specific pattern for using the wells on the plate. A proposed layout of standards (Std) and samples (S) to be measured in duplicate is shown below.

NOTE: If sample blanks are included in the assay, it is necessary to reserve some wells of the plate for these blanks

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 1	Std 1	S1	S1	S9	S9	S17	S17	S25	S25	S33	S33
B	Std 2	Std 2	S2	S2	S10	S10	S18	S18	S26	S26	S34	S34
C	Std 3	Std 3	S3	S3	S11	S11	S19	S19	S27	S27	S35	S35
D	Std 4	Std 4	S4	S4	S12	S12	S20	S20	S28	S28	S36	S36
E	Std 5	Std 5	S5	S5	S13	S13	S21	S21	S29	S29	S37	S37
F	Std 6	Std 6	S6	S6	S14	S14	S22	S22	S30	S30	S38	S38
G	Std 7	Std 7	S7	S7	S15	S15	S23	S23	S31	S31	S39	S39
H	Std 8	Std 8	S8	S8	S16	S16	S24	S24	S32	S32	S40	S40

Example of plate layout for the Cysteine Assay Kit

7. Sample preparation

The following sample preparation protocols are intended as a guide only. The optimal conditions for sample preparation must be determined by the end user. It is recommended to use fresh samples. If it is not possible, aliquot and store samples with minimal freeze/thawing.

BQC Cysteine Assay Kit can be used to determine Cys content in a wide variety of samples like animal and plant tissues, cells, plasma, serum and other biological fluids.

Animal Tissues. Weigh 0.1 g tissue, add 1 mL of **Reagent D** and homogenize on ice. Centrifuge at 11000 g for 10 minutes at 4 °C. Use the supernatant and place it on ice to be tested.

Plant Tissues. Weigh 0.1 g tissue. Add 1 mL of **Reagent D** and mash. Ultrasonic break in ice bath for 5 minutes and centrifuge at 11000 g for 10 minutes at 4 °C. Use the supernatant and place it on ice to be tested.

Cells and bacteria. Collect $5 \cdot 10^6$ cells into the centrifuge tube, wash cells or bacteria with cold PBS and discard supernatant after centrifugation. Add 1 mL of **Reagent D** to ultrasonically disrupt cells or bacteria for 5 minutes. Centrifuge at 11000 g for 10 minutes at 4 °C. Use supernatant for the assay and place it on ice to be tested.



Serum, plasma, and other liquid samples. Dilute the sample with **Reagent D** to fit within the calibration curve. If it is necessary, centrifuge at 11000 g for 10 minutes at 4 °C and use supernatant for assay. Keep on ice to be tested.



Reagents and materials required for sample preparation are not supplied with the kit. Before doing sample preparation, consider the volume of sample required per test; see **Technical specifications** section.



Make sure that interfering substances present in the sample do not give a significant background. Run proper blanks as necessary. It is recommended to assay different sample dilutions to ensure the values fall within the standard curve.

8. Assay protocol

Prepare and mix all reagents thoroughly before use. Each standard, sample or blank should be assayed at least in duplicate.

- 1  Set up the plate design
- 2  **Standard** and **sample** wells: add **20 µL** of **sample** or **standard** to each well

! **Optional:** Sample blank wells: Add **20 µL** of **sample**
- 3  Add **100 µL** of **Reagent A** to each well
- 4  Add **100 µL** of **R.B Working Solution** to each well


! **Optional:** Sample blank wells: Add **100 µL** of **ddH₂O** or **buffer**
- 5  Kept at **RT** for **15 minutes**
- 6  Read the **absorbance** of all wells at **600 nm**

If you need to **adapt this kit** for another form of the assay (for example cuvette), **contact us at** info@bioquochem.com

9. Data analysis

ANALYSIS OF THE STANDARDS

- Calculate the average absorbance of the standards.
- Subtract the average absorbance of the reagent blank (**Std 8**) from the average absorbance of all the standards to obtain the blank-corrected absorbance of the standards.
- Create a standard curve by plotting the blank-corrected absorbance of the standards as a function of the standard concentration (see **Standard Preparation** section). A typical standard curve ($y = \text{slope} \cdot x \pm \text{intercept}$) for this assay is shown below.



Standard curve for Cysteine Assay Kit

- ! This standard curve is an example of the data typically obtained with this kit. **DO NOT USE** this standard curve to calculate the Cys content of your samples. A new standard curve must be performed by the end user.

ANALYSIS OF THE SAMPLES

- Calculate the average absorbance of the samples.
- Subtract the average absorbance of the reagent blank (**Std 8**) from the average absorbance of each sample to obtain the blank-corrected absorbance of the samples (**A_s**).
 - ❗ If sample blanks are assayed and they are significant, subtract them from the blank-corrected absorbance of the samples
- Calculate the Cys content from a sample as mM Cys using the equation obtained from the linear regression of the standard curve by substituting the blank-corrected absorbance for each sample (**A_s**).

$$\text{Cysteine (mM)} = \left(\frac{A_s - \text{intercept}}{\text{slope}} \right)$$

When working with diluted samples the concentration values obtained must be multiplied by the dilution factor to obtain the Cys (mM) value of the undiluted sample.

10. Troubleshooting

This troubleshooting table provides potential sources and solutions for common problems observed with BQC Assay Kits. **The problems listed below could occur when using any BQC Assay Kit.** They are not specific for this Assay Kit.

Problem	Possible Cause	Recommended Solution
Wells have color but there is no reading	Plate read at incorrect wavelength	Check the wavelength used in the assay
	Incorrect microplate	Use the correct microplate for your application UV/Vis: transparent Fluorescence: black wells/transparent bottom
Standard readings do not follow a linear pattern	Pipetting errors in preparation of standards	Avoid pipetting small volumes (<5 µL) Be careful not to splash from well to well
	Air bubbles formed in well(s)	Use reverse pipetting technique
	Standard stock is at incorrect concentration	Always refer to dilutions described in Assay preparation
	Improperly thawed reagents	Thaw all components completely and mix well before use
	Use of improperly stored reagents	Store the components appropriately Use fresh components from the standard curve
	Incorrect incubation times or temperatures	Refer to Refer to Assay protocol
Dispersion of standard and sample readings	Pipetting errors	Avoid pipetting small volumes (<5 µL) Be careful not to splash from well to well
	Air bubbles formed in well(s)	Use reverse pipetting technique

Problem	Possible Cause	Recommended Solution
Sample erratic values	Samples contain interfering substances	Dilute sample further (if possible)
	Inappropriately stored samples or samples used after multiple freeze-thaw cycles	Use fresh samples or store appropriately until use
	Samples not deproteinized	Use an appropriate deproteinization protocol
	Cells/Tissue samples not homogenized completely	Repeat the sample homogenization
	Inappropriate sample dilution buffer	Refer to Assay preparation
Sample reading fall outside the detection range	Samples are too diluted/concentrated No analyte/activity is observed in the sample	Re-assay using different sample dilutions

STILL HAVING PROBLEMS?

Contact BQC if you have any further questions, our team will be pleased to help you:



Phone + 34 985 26 92 92



E-mail info@bioquochem.com



Business hours Monday-Thursday: 8.30 to 17.00 (CEST)
Friday: 8.00 to 15.00 (CEST)

11. Additional information

Cysteine Assay Kit is a simple and reproducible (RSD< 4%) assay for determining Cys content in a wide variety of biological samples.

To calculate Cys content by protein concentration, **BQC BCA Assay (KB03005)** or another Protein Quantification Assay Kit can be used to determine the total amount of protein in the sample.

If unexpected results are obtained running your samples, please contact us at info@bioquochem.com

12. Related products

More products available at bioquochem.com

Reference	Product
KB03005	BCA Protein Quantification Assay Kit
KB03031	Protein Concentration Assay Kit
KB03008	Protein Carbonyl Assay Kit

13. Warranties and limitation of liability

BQC shall not in any event be liable for incidental, consequential or special damages of any kind resulting from any use or failure of the products, even if BQC has been advised of the possibility of such damage including, without limitation, liability for loss of use, loss of work in progress, downtime, loss of revenue or profits, failure to realize savings, loss of products of buyer or other use or any liability of buyer to a third party on account of such loss, or for any labor or any other expense, damage or loss occasioned by such product including personal injury or property damage is caused by BQC's gross negligence. Any and all liability of BQC hereunder shall be limited to the amounts paid by the buyer for the product.

Buyer's exclusive remedy and BQC's sole liability hereunder shall be limited to a refund of the purchase price, or the replacement of all material that does not meet our specifications.

Said refund or replacement is conditioned on buyer giving written notice to BQC within 30 days of shipment.

Expiration date: 1 year from the date of fabrication. Expiration date is indicated on the outside of the box.

For further details, please refer to our website [bioquochem.com](https://www.bioquochem.com)



Edificio CEEI | Parque Tecnológico de Asturias,

33428 Llanera, Asturias

Info@bioquochem.com

[linkedin.com/bioquochem](https://www.linkedin.com/company/bioquochem)

www.bioquochem.com