# **GeneJect**<sup>™</sup>

## Vesatile transfection reagent

www.aabiot.com/en/geneject



## **Advantages**

- Suitable for DNA, siRNA and miRNA transfection
- Low cell toxicity
- Efficient transfection of wide variety of primary cells and cell lines

## Description

GeneJect<sup>™</sup> is a transfection reagent, which is suitable for the delivery of plasmid DNA and short oligonucleotides, such as siRNAs as well as miRNA mimics and inhibitors.

GeneJect  $^{\rm M}$  is less toxic when compared with lipofection. Toxicity is higher at low cell confluency.

GeneJect<sup>™</sup> is efficient in the transfection of common cell lines of different origin, including HeLa, HEK293, PPC-1, HaCaT, A549, BEAS-2B, CCD19Lu, THP-1 and KOP (fetal bovine oropharynx cells). It also enables transfection of low level transfection cells including MDBK (Madin Darby Bovine Kidney cells).

GeneJect<sup>™</sup> is efficient in the transfection of wide variety of primary cells, such as human primary keratinocytes, primary skin fibroblasts, melanocytes and monocytes.

100  $\mu l$  of GeneJect^\* allows for up to 50 transfections in 12-well tissue culture plates format.

#### Uptake of GeneJect<sup>™</sup> - plasmid DNA nanocomplexes

Cells were transfected with pEGFP plasmid DNA on a 12 well plate with 1 µg DNA and 2 µl GeneJect<sup>™</sup> per well and visualised under fluorescence microscope after 24h from transfection.



**HEK** - Easy to transfect Human embryonic kidney cells. **MDBK** - Hard to transfect Madin Darby bovine kidney cells.

#### Uptake of GeneJect<sup>™</sup> - miRNA nanocomplexes

Dy547-labelled (red) miRNA mimics were transfected with either GeneJect<sup>™</sup> or with a lipofectamine type of reagent siPORT-NeoFX at 30 nM concentration and stained after 24 h with DAPI (blue) and Alexa Fluor A488 Phalloidin (green).



Dy547 Phalloidin DAPI

**A549** - Adenocarcinomic human alveolar basal epithelial cells, **BEAS-2B** - Human bronchial epithelial cell line, **CCD19-Lu** - Human lung fibroblast cells, **PPC-1** - human prostatic cell line.

### siRNA Transfection performance of GeneJect™

#### Inhibition of target gene (IRAK1) by GeneJect<sup>™</sup> delivered miR-146a mimic

Cells were transfected either with control (cont) or miR-146a mimics at 30 nM concentration.

After 48 hours, a total RNA was purified and analyzed by RT-qPCR using primers specific to IRAK1. Relative expression compared to cont in each cell type (=1) is shown. Data represent mean  $\pm$  SEM of four independent transfection. Student's t-test, \*\* P < 0.01.



**PPC-1** - Human prostatic cell line, **A549** - Adenocarcinomic human alveolar basal epithelial cells, **KC** - Human primary keratinocytes, **THP-1** - Human monocytic cell-line.

## GeneJect<sup>™</sup> efficiency in siRNA delivery

Human primary keratinocytes were transfected either with control (cont) or indicated specific siRNAs at 30 nM concentration using GeneJect<sup>™</sup>. After 48 hours, a total RNA was purified and analyzed by RT-qPCR using primers specific to IRAK1 and CARD10.

Relative expression compared to cont siRNA (=1) is shown. Data represent mean  $\pm$  SEM of three independent transfection. Student's t-test, \* P < 0.05. \*\* P < 0.01.

