

# TCA DEPROTEINIZING ASSAY KIT

## KB03026-100/200/400 Tests

#### **DESCRIPTION AND USE**

Proteins may interfere with some assays, affecting accuracy and sensitivity. When ultrafiltration cannot be done, other chemical removal alternatives can be considered. BQC TCA Deproteinizing Kit is recommended for the deproteinization of samples prior to assaying low molecular weight metabolites. The precipitated proteins will remain nonfunctional.

The BQC TCA Deproteinizing Kit ensures a protein removal efficiency over 99 % with very low sample dilution that includes a neutralizing solution to adjust the pH. The volume of sample required per test is 150 µL (samples with low protein concentration) or 90 µL (samples with high protein concentration)

#### **MATERIALS SUPPLIED**

Item	No. Tests	Quantity
TCA Solution	100	1
	200	2
	400	4
Neutralizing Solution	100	1
	200	2
	400	4

#### STORAGE AND STABILITY

On receipt store kit components at RT. Do not use after the expiration date stated on the packaging.

#### RELATED PRODUCTS

Product	Reference
Bradford Protein Assay Kit	KB03003
ABTS Antioxidant	KF01002
Capacity Assay Kit	

### **ASSAY PROTOCOL**

1	10 min	Place the solutions on ice to ensure they are cold	
2		In a microtube, <b>mix</b> your <b>sample</b> with the <b>TCA Solution</b> For samples with <b>high protein concentration</b> use a <u>6:1 ratio</u> . For example: 90 µL of sample with 15 µL of TCA Solution. For samples with <b>low protein concentration</b> use a <u>10:1 ratio</u> . For example: 150 µL of sample with 15 µL of TCA Solution.	
3	1 min	Vortex	
4	15 min	Keep microtubes on ice	
5	10 min	Centrifuge at 10000 x g at 4 °C	
6		Collect the supernatant in other microtube. If proteins are required, collect the pellet, and freeze at -80 °C	
7		Add 10 µL of the Neutralizing Solution.	
8		Check that the pH is neutral with a pH paper test. If necessary, adjust to pH 7 with the Neutralizing Solution.	
9		Assay directly or freeze at -80 °C until the day of the assay	

For future experiments and calculations consider that the sample is diluted throughout the deproteinizing assay protocol. Consider the dilution factor performed when analyzing the results.

FOR RESEARCH USE ONLY

Explore our web











