

ZINC SULPHATE DEPROTEINIZING ASSAY KIT

KB03028-10/25/50 Tests

DESCRIPTION AND USE

Proteins may interfere with some assays, affecting accuracy and sensitivity. When ultrafiltration cannot be done, other chemical removal alternatives can be considered. BQC Zinc Sulphate Deproteinizing Kit is recommended for the deproteinization of samples prior to assaying low molecular weight metabolites.

The BQC Zinc Sulphate Deproteinizing Kit ensures a protein removal efficiency over 99 % with no sample dilution and do not require pH correction. The volume of sample required per test is 400 µL. The number of vials provided equals the number of samples to be used.

MATERIALS SUPPLIED

Item	No. Tests	Quantity
Zinc Sulphate	10	10
	25	25
	50	50

STORAGE AND STABILITY

On receipt store kit components at RT. Do not use after the expiration date stated on the packaging.

RELATED PRODUCTS

Product	Reference
Catalase Activity Assay Kit	KB03012
MDA-TBARs Assay Kit	KB03016

ASSAY PROTOCOL

1	10 min	Place the solutions on ice to ensure they are cold	
2		Add 400 μL of your sample to a vial. Use 1 vial per sample.	
3	1 min	Vortex	
4	15 min	Keep microtubes on ice	
5	10 min	Centrifuge at 10000 x g at 4 °C	
6		Collect the supernatant in other microtube. If proteins are required, collect the pellet, and freeze at -80 °C	
7		Assay directly or freeze at -80 °C until the day of the assay	

FOR RESEARCH USE ONLY

Explore our web www.bioquochem.com











