## Ready-Strain Cell Straining Kits



- Includes strainer, 50mL tube, & cap
- Opens directly under the hood reduces contamination risks
- Preassembled ready to use
- 💍 Guaranteed sterile gamma
- Packaged in easy-open "tear-strip" bags

ReadyStrain<sup>™</sup> is the first product of its kind, providing all the necessary parts for filtration of cell suspensions and combining them into a **single**, **preassembled**, **ready to use**, **sterile packet**. These kits eliminate the contamination risks associated with opening multiple packages, and assembling separate pieces (i.e. tubes strainers and caps). The individual sterile packets mean researchers are no longer forced to expose entire bags of 50mL tubes and corresponding caps in order to perform one filtration. Instead they can open a single ReadyStrain packet directly under the hood.

Each kit contains a *sterile* SureStrain™ Strainer, pre-attached to a *sterile* 50mL tube, with a separate *sterile* screw-cap ready for assembly, processing and storage after the filtration is complete. Kits are supplied in easy open "tear-strip" bags and gamma sterilized.



2600A Main Street Extension Sayrevill, NJ 08872 Ph: 908-279-0280 www.mtcbio.com

## What's included:

1 x 50mL Centrifuge Tube w/ graduations

1 x Attached SureStrain Cell Strainer (40, 70 or 100µm)

1 x Screw Cap for 50mL Centrifuge Tube

Kits supplied sterile, preassembled, in individual easy-open bags, 50/pk



## Also available: PolyPestle<sup>™</sup> Dual-sided Pestle



Cat. No.	<u>Description</u>	<u>Qty</u>	<u>Price</u>
C5040	Cell straining kit w/ strainer, 50mL tube and screw-cap, 40µm, blue, sterile, individually wrapped	50	\$139
C5070	Cell straining kit w/ strainer, 50mL tube and screw-cap, 70µm, white, sterile, individually wrapped	50	\$139
C5100	Cell straining kit w/ strainer, 50mL tube and screw-cap, 100µm, yellow, sterile, individually wrapped	50	\$139
C4010	PolyPestle <sup>™</sup> dual-sided pestle for strainers & 1.5mL microtubes, sterile, individually wrapped	50	\$39