



Transport Swab


Environmental Sample
Collection and Dilution Device

Hygiena™ Transport Swab is the ideal sampling device for customers looking for an environmental sample collection device containing more than 1mL of enrichment media. Transport Swab makes sample collection efficient and easy.

Benefits:

- Ready-to-use and convenient—pre-filled bottles come sterilized and have ample head space for additional sample
- All-in-one sample collection system
- More broth for high-volume samples
- Available in 4 mL and 10 mL formats
- Available in Lethen or Buffered Peptone Water
- Perfect for plates or dried media film
- Collect your samples by swabbing the surface and recapping
- 18-month shelf life at refrigerated temperature (2 – 8 °C)

Transport Swab Procedure



Catalog No.	Description	Quantity
HTS-4	Transport Swab (Lethen 4 mL)	100
HTS-10	Transport Swab (Lethen 10 mL)	100
HTS-10BPW	Transport Swab (Buffered Peptone Water 10 mL)	100



Find support documents, instructional videos, and more at www.hygiena.com