

Filtration Accessories



PES / PC Vacuum Bottle

Specification	1200 ml, PES vacuum bottle	1200 ml, PC vacuum bottle
Catalog number	197000-11	197000-11-PC
Graduated increment	100 ml	100 ml
Overflow protection	Yes	Yes
Cap with suction fitting	Yes	Yes
Hose barb	ID8	ID8
Autoclavable	Yes	Yes



PP Heavy-Duty Vacuum Bottle **NEW**

Specification	4000 ml, PP vacuum bottle	4000 ml, PP vacuum bottle with drain
Catalog number	197200-47	197200-48
Connector	Hose barb	Hose barb
Graduated increment	500 ml	500 ml
Overflow protection	Yes	Yes
Hose barb	ID4 ~ ID8	ID4 ~ ID8
Drain device	-	Yes
Autoclavable	Yes	Yes



Multilayer Filtration Adaptor **NEW**

Specification	Multilayer filtration adaptor
Catalog number	167230-05
Main material	Borosilicate glass
Material of membrane support	SS screen
Filter diameter	47 mm
Filter holder	300 ml glass filter holder (optional)



Continuous Filtration Kit **NEW**

Specification	Continuous filtration kit
Catalog number	167210-36
Main material	PTFE
Material of hose barb	PTFE
Filter holder	300 ml glass filter holder (optional)
Hose barb	ID6



* PES : Polyethersulfone, PC : Polycarbonate, PP : Polypropylene

Filtration Accessories



Membrane Filters

Specification	Glass fiber filter	MCE membrane
Catalog number	167100-20	167100-52
Brand	Pall	Pall
Filter diameter	47 mm	47 mm
Pore size	1µm	0.45µm
Material of filter	Glass fiber	MCE
Package	100/PK	200/PK



GL45 Suction Cap / Filtration Adaptor

Specification	GL45 suction cap	GL45 filtration adaptor
Catalog number	197000-60	197000-65
Material of cap	PP	PP
Material of gasket	-	PTFE
Stopper	-	No.8 stopper
Hose barb	ID6	ID4 ~ ID8
Autoclavable	Yes	Yes

Application for GL45 Suction Cap

Waste Bottle

Suitable for collecting waste media by using GL45 bottle in connection with suction pump and accessories.




Moisture Trap Bottle

Available to absorb moisture with desiccant which turns pink when saturated.



Vapor Absorption (Neutralization) Bottle

Ideal for absorbing acid or alkaline vapor to prevent harmful vapors from spreading into the environment or pump.



Safety Bottle

Works as a safety bottle to prevent liquid from being sucked into pump.



* Please refer to P.66 to select GL45 Laboratory Glass Bottle