

Aluminum Multi-Branch Manifold



VIDEO

2
YEAR
WARRANTY

Features

▶ Light and handy aluminum design

All parts of MultiVac 310 / 610 manifold are made of aluminum which is light, cost-effective and autoclavable.

▶ Quick-Fit hose barb

Detachable hose barb can be connected to either side of the manifold quickly towards the vacuum source.

▶ Individual control valve

Each branch permits individual control.

▶ Compatible with different sizes of filter holders

MultiVac 310 / 610 - MS manifold includes No.8 silicone stoppers which can accept various types of filter holders at option.

▶ Two-Year Warranty

Applications

- ▶ Suspended solid test
- ▶ Microbiological test
- ▶ Vacuum filtration

Ordering Information

180310-01

MultiVac 310 - MS, 3 Places Aluminum Manifold with No. 8 stopper x 3

180610-01

MultiVac 610 - MS, 6 Places Aluminum Manifold with No. 8 stopper x 6


200300-01

MF 3, 300 ml Magnetic Filter Holder with lid kit

167110-16

No.8 Silicone Stopper

MultiVac 310 - MS



Quick-Fit hose barb

MultiVac 610 - MS



Specification

Model	MultiVac 310 - MS	MultiVac 610 - MS
Number of branch	3	6
Filter diameter	Optional	Optional
Funnel capacity	Optional	Optional
Effective filtration area	Optional	Optional
Hose barb	5/16 inch (8 mm)	5/16 inch (8 mm)
Dimension (LxWxH)	48 x 13 x 16 cm	85x 13 x 16 cm
Part	Material	Material
Body	Aluminum	Aluminum
Control valve	Aluminum	Aluminum
Handle	Aluminum	Aluminum
Hose barb	Aluminum	Aluminum